## Goldmann AT 900 / 870

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## Introduction

We would like to thank you for your decision to purchase this Haag-Streit product.

If the instructions in this manual are carefully followed we are confident that this product will give you reliable and trouble-free usage.

## Purpose of use

The **Applanation Tonometer 900** is an accessory to the Slit Lamp 900® which serves to measure ocular pressure and can be fitted to Haag-Streit slit lamps.

Due to its versatility, the Applanation Tonometer can also be mounted onto other slit lamps.

The **Applanation Tonometer 870** can also be used as an accessory for the slit lamps of other manufacturers.

**Applanation Tonometers** are used, at room temperature, in the examination, diagnosis and documentation of the human eye.

Applanation tonometry may only be performed by healthcare professionals who are qualified by their respective authorities to do so.

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• The Goldmann tonometer is still the gold standard. But this does not prevent the user from studying new publications on the influence of biomechanical parameters of the cornea on the measurement of the intraocular pressure.

## Contents

<b>1</b>	<b>Safety</b> .....	5
<b>2</b>	<b>Overview</b>	
2.1	General structure .....	11
2.2	Which slit lamp model with which tonometer ..	11
2.3	Description of models .....	12
2.4	Accessories .....	12
2.5	How the pression is measured .....	15
<b>3</b>	<b>Operating</b>	
3.1	Preparing the patient .....	17
3.2	Preparation of the equipment .....	17
3.3	Instructions to the patient .....	19
3.4	Correct measurement .....	21
3.5	Sources of error .....	23
3.6	Tips for the measurement .....	31
3.7	Astigmatism .....	31
<b>A</b>	<b>Installation</b>	
A.2	Tonometer AT 900 M/Q .....	33
A.1	Tonometer AT 900 C/M .....	33
A.3	Tonometer AT 900 BQ .....	33
A.4	Tonometer AT 870 .....	33
<b>5</b>	<b>Maintenance</b>	
5.1	Check tonometers AT 900 .....	35
5.2	Check tonometer AT 870 .....	39
5.3	Cleaning / disinfecting of measuring prisms ..	43
5.4	Accessories disinfecting / storing .....	43
<b>6</b>	<b>Technical specifications</b> .....	45
Short	<b>How to disinfect?</b> Maintain the tonometer ...	48

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# 1 Safety

## Ambient conditions

- Transportation Temperature -40 °C to +70 °C  
Air pressure 500 hPa to 1060 hPa  
Relative humidity 10% to 95%
- Storage Temperature -10 °C to +55 °C  
Air pressure 700 hPa to 1060 hPa  
Relative humidity 10% to 95%
- Working Temperature +10 °C to +35 °C  
Air pressure 800 hPa to 1060 hPa  
Relative humidity 30% to 75%

## Installation, dispatch

- To avoid condensation, allow the instrument to adjust to room temperature for several hours before unpacking.
- The Applanation Tonometer should be verified for damage after unpacking.
- Installation is to be carried out according to Appendix A and functional verifications according to Appendix B.
- Check: do the connection pieces fit (tonometer on slit lamp, measuring prism)?
- Defective equipment should always be returned in an appropriate packaging.

## Operation, surrounding


- Only qualified and trained personnel should operate the equipment, the training is at the owner's responsibility.
- Use only Haag-Streit accessories.
- Avoid humidity.

## Cleaning of housing

- Use a slightly moistened cloth.
- Soap can be employed.
- No liquids or corrosive agents.
- Avoid direct spraying on the equipment.

## Warranty, product liability

- The instrument should be operated in accordance with the chapter 'Safety'. Incorrect operation can damage the instrument. Thus no warranty claims can be accepted.
- In particular, it is forbidden to loosen any screws; this may reduce the accuracy of the measurement.
- Continued use of an instrument that has been damaged by incorrect operation can lead to personal injury. The manufacturer cannot accept liability in this case.
- Repairs and alterations on this equipment should only be carried out by Haag-Streit. Contact your Haag-Streit distributor.

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
**Avoid examinations in case of eye infections or injured corneas**

- There is no contraindication for performing tonometry. Medical assessment is necessary and due care is to be exercised.


**Only clean and disinfected measuring prisms should be used**

- Examinations should only be conducted with cleaned and disinfected measuring prisms.
- Tonometer prisms are not dispatched in a disinfected state and should always be disinfected before use.
- The disinfecting of contact glasses eg. with cotton wool or damp cloth without putting them into disinfecting fluid is insufficient and not allowed.
- Incorrect medical disinfecting fluid can cause patients to be infected and can also cause damage to the eyes.
- Residue of the disinfecting may cause provoking or caustic reactions on the patient's eye, therefore rinse thoroughly.
- Measuring prisms must be removed from the holder and be put into disinfecting fluid.
- Cleaning and disinfecting is to be carried out according to Appendix B, Paragraph B.3 and short manual.

**No damaged prisms should be used**



- **Before each use** the contact surface of all prisms should be verified for contamination or damage (such as scratching, rips or sharp edges). The slit lamp microscope is best used for this purpose at a magnification of 10x to 16x.
- If cracks form, the disinfectant may get into the hollow interior and can enter the patient's eye while tonometry is taking place, leading to irritation of the cornea. Defective measuring prisms should not be used, in order to avoid any damage to the cornea of the patient.
- In case of inaccurate measuring results, please observe the functional verifications according to Appendix B.

**Control**

- This type of sensitive measuring equipment is to be checked after exposure to external force (e.g. unintentional shocks or being allowed to fall) in accordance with Appendix. If repair is necessary please contact your Haag-Streit distributor.

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
**Measuring prism - what time of usage?**

- Due to the number of variables involved (type and concentration of disinfecting material used, number of patients, handling, etc.) it is almost impossible to put a figure on the number of times and / or period of time that a measuring prism can safely be used.
- Haag-Streit recommend that no measuring prism may be used longer than two years under 'normal' conditions, i.e. as per the instructions contained within this instruction manual.
- A measuring prism may be retained indefinitely as a spare part under room temperatures and protected from sun exposure, deterioration does not start until it is used for the first time.
- This recommended time is overruled if the prism is in any way damaged within this period of time; damaged prisms should be replaced immediately.

**Statutory requirements**


- The Applanation Tonometer is designated as Class I product with measuring functions under the EEC Directive 93/42 for medical equipment products.
- The CE designation is being then granted on the basis of this Directive as per Conformity Module B.
- The Applanation Tonometer complies with the Tonometer regulation ISO 8612.
- A copy of the declaration of conformity of the present instrument can be requested at any time from Haag-Streit.
- Statutory accident regulations are also to be observed.
- National calibrating regulations for measuring equipment are also to be observed.

• Classification

CE-Regulation 93/42 EEC    Class Im (measuring functions)  
 FDA                                    Class II

**Pictograms**

- a Read very carefully (in the instruction manual)!
- b Legend plate

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


## 2 Overview


### 2.1 General structure

- 1 Measuring prism
- 2 Feeler arm
- 3 Legend plate (underside)
- 4 Control weight insert
- 5 Revolving knob with measuring drum


### 2.2 Which slit lamp with which tonometer?

Slit lamp	Tonometer model AT			
	900 C/M	900 M/Q	900 BQ	870
• Haag-Streit				
BC 900®	●			
BD 900®	●			
BM 900®	●	●		
BM 900® V	●	●		
BP 900®		●	●	
BQ 900®		●	●	
BX 900®		●	●	
• Others				●

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### 2.3 Description of models

The **Applanation Tonometer AT 900 C/M (6)** is produced for those who wish the tonometer to remain permanently on the slit lamp. It is mounted on a pivot on the microscope and for examination is swung forward in front of the microscope. The applanation surface is to be examined in a monocular manner through the left ocular.



The **Applanation Tonometer AT 900 M/Q (7)** is mounted on the guide plate over the slit lamp axis during tonometric examinations. For examinations through the right or left ocular, the appliance can be mounted in two positions on the guide plate. The feeler arm protrudes from below with prism measuring into the beam path of the microscope and the illumination unit.

Apart from the working position, the **Applanation Tonometer AT 900 BQ (8)** has also two resting positions.

For perfect tonometric examination the tonometer is swung in front of the microscope and the illumination unit moved from the left to touch the bracket of the tonometer. In this position, the left and right eye of the patient can be examined without difficulties (no 60° position).


The applanation surface is to be examined in a monocular manner through the right ocular of the stereo-microscope.

The **Applanation Tonometer AT 870** is to be mounted over the microscope of the slit lamp. The feeler arm protrudes from above into the path of the microscope and the illumination unit. The tonometer is ready for use with a single swing movement of the feeler arm.

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## 2.5 How the pressure is measured

*The Applanation Tonometer functions according to the 'Goldmann method': The measuring of the pressure required to maintain a uniform applanation of the surface of the cornea.*

The cornea is applanated with a measuring prism of organic glass which is fitted in a ring shaped holder at the end of the feeler arm. The contact surface of the measuring prism has a diameter of 7.0 mm. It is planar with rounded edges to obviate any damage to the cornea.

The measuring prism is brought into contact with the eye of the patient by moving the slit lamp forwards. The measuring drum is turned to increase pressure on the eye until a continuously regular surface of 3.06 mm  $\phi$  = 7.354 mm<sup>2</sup> has been applanated.

Position of the measuring drum	Force	Pressure
1	9.81 mN	1.33 kPa $\approx$ 10 mm Hg

The intraocular pressure in mm Hg is found by multiplying the drum reading by ten.

According to ISO 1000 the derived unit of force is the N (Newton) and its decimal sub-multiple the mN (milli Newton).

The conversion of mm Hg into the SI pressure unit kPa (Kilo Pascal) can be made on the *comparison scale* (10).

The measurement of the flattened surface is made directly on the cornea. The built-in duplication system in the measuring prism divides the image and displaces the two semi circular halves from each other by 3.06 mm.


Before measuring, it will be necessary to anaesthetise the cornea locally, to place a strip of fluorescein paper in the conjunctival sac and switch on the blue filter of the slit lamp. The inner border of the ring represents the line of demarcation between the cornea flattened by applanation and the cornea not flattened.

The major advantage of applanation tonometry is the small bulbus deformation which amounts to only 0.56 mm<sup>3</sup>. The values found by this method of tonometry are only slightly influenced by scleral-rigidity and radius of corneal curvature, the intraocular tension increases by about 2.5% only.

The principle of applanation tonometry is simple. The careful construction of the apparatus is a guarantee of its continued correct performance.

Nevertheless, in order to obtain exact results, the method of use indicated must be strictly followed.

The readings obtained are based on the assumption that a cornea with a 'normal' thickness is being measured, alteration of the corneal thickness leads to changes in the measured IOP. A 'normal' corneal thickness is considered to be within the area of 530 to 560 microns.

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### Magic formula?

*There is no universally agreed formula for the adjustment of the measured tonometry reading against altered corneal thickness.*

*The users should bear these two facts in mind and keep themselves informed, through the current literature, with regard to the various formulae for adjusting measured tonometry readings in the event of thin or thick corneas being measured.*

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### 3 Operating

#### CAUTION

It is imperative to read the chapter 'Safety' and to observe its precautions before operating the equipment.

#### NOTICE

Perfect functioning can only be guaranteed with the use of an original Haag-Streit prism.

Sterilised fluorescein paper strips should always be used, as pathogenic excitors thrive well in fluorescein solutions.

Should the measuring prism come into contact with the cornea without the drum having previously been correctly set, vibrations can occur in the feeler arm, which will produce unpleasant feelings for the patient.

#### 3.1 Preparing the patient

1. Both eyes must always be anaesthetised (e.g. 2-3 drops each of an anesthetic within half a minute) to reduce movements of the lids during examination.

2. Place a fluorescein paper strip near the external canthus in the lower conjunctival sac. After a few seconds the lacrimal fluid is sufficiently coloured and the paper can be removed.

When using drops, a solution of sodium fluorescein 0.25 % to 0.5 % is advisable. Should you wish to instil a solution of 1 % or 2 %, introduce a small amount of liquid into the conjunctival sac by means of a glass rod.

3. Place the head of the patient on the chin rest.

#### 3.2 Preparation of the equipment

##### For all Haag-Streit slit lamps and tonometers

4. Before measuring, make sure that the eyepieces are correctly focused.

5. Set the magnification at 10x.


6. Set the power supply / the potentiometer of the illumination to minimum position.

7. Bring the blue filter into the beam of the slit lamp and open the slit diaphragm completely.


8. The feeler arm should be so engaged that the axis of the measuring prism and the microscope are convergent.

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

## AT 900 C/M and AT 900 MQ

10. Swing the illumination unit to the left.


11. Bring the tonometer from the right in front of the microscope into the measuring position located by a notch.

12. Swing the illumination unit from the left to touch the holder of the tonometer. This is the only position of the illumination unit permitting perfect tonometric examination of the left and right eye of the patient (no 60° position). This position facilitates the opening of the lids of the patient by the examiner where this procedure is necessary for measurement. The illumination of the applanated surface is practically free of reflections through the measuring prism.

Examination

- with AT 900 C/M through the left ocular
- with AT 900 BQ through the right ocular


AT 900 C/M

AT 900 M/Q

AT 900 BQ

AT 870

## AT 900 M/Q


10. During examinations with the tonometer, the angle between the illumination unit and the microscope in the left or right ocular, should be 60° approximately, so that the image appears bright and free of reflections. Other solution: illumination from the back at about 10°.

AT 900 C/M

AT 900 M/Q

AT 900 BQ

AT 870

## AT 870

10. Swing the pressure arm with the prism attached into the beam path of the illumination unit and the microscope.

11. Bring the blue filter into the beam path of the illumination unit in the slit lamp and open the slit diaphragm completely. The angle between illumination unit and microscope should be 60° approx., so that the image appears bright and free of reflections. Other solution: illumination from the back at about 10°.

12. The microscope is to be adjusted before the examination and according to the refraction divergence of the examiner, so that the fluorescein rings are observed sharply focused during the examination.

13. Turn the measuring drum to setting 1.


14. Use average illumination strength.

AT 900 C/M

AT 900 M/Q

AT 900 BQ

AT 870

### 3.3 Instructions to the patient

1. The face of the patient should be firmly held against the chin and forehead rests. A headband to hold the head can also be used.

2. Instruct the patient to look straight ahead. If necessary use the fixation target to steady the eyes.

3. It is recommended that the patient is repeatedly asked to keep his eyes wide open during the examination. If need be the examiner may hold open the lids of the examined eye with his fingers provided that no pressure is applied to the eye.

The angle between the microscope and the illumination unit is to be reduced to around 10° also on models AT 900 M/Q and AT 870, whereby the light beam will go through the body of the prism.

A reflection-free image should be ensured.

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### 3.4 Correct measurement

1. Immediately before taking the measurements, the patient should be made to close the eyes briefly so that the cornea becomes sufficiently moistened with the lacrimal fluid and fluorescein.

2. By moving the slit lamp forward, the measuring prism is brought into contact with the centre of the cornea over the area above the pupils.

The limbus of the cornea will be illuminated with a bluish hue. This illumination is best observed by direct sight from the opposite side of the illumination unit by the examiner.

As soon as the limbus of the cornea illuminates, cease any further forward movement of the slit lamp immediately.

3. After obtaining contact, commence observing the cornea through the microscope. A regular pulsation of the two semi circular fluorescein rings, which can be of varying size dependent on the ocular pressure, when the measuring drum is at setting 1, will show that the tonometer is in the correct measuring position.

The requisite correction is made with the guide lever until the applanated surface is observed as two semi circular surfaces of the same size in the centre of the field of sight (11).

Small lowering adjustments of the slit lamp with the guide lever have no effect on the size of the semi circular images.

4. The pressure on the eye is increased by turning the measuring drum on the tonometer, until the edges of both fluorescein rings just meet = **correct position** (12).

The edges cross over each other with the pulsation of the eye.

The width of the fluorescein ring around the contact position of the measuring prism should be 1/10 approximately of the diameter of the applanation surface (0.3 mm).

5. Reading the scale:


- The reading
- multiplied by 10
- gives the ocular pressure in mm Hg

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
*Notice* The following are schematically simplified illustrations!

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### 3.5 Sources of error

#### Images in the ocular


Wrong fluorescein ring . . . . .	13 - 14
Wrong distance to the patient. . . . .	15 - 16
Position too much at right / left . . . . .	17 - 20
Position too high / low. . . . .	22 - 25
Inaccurate pressure . . . . .	27 - 29

#### The fluorescein ring is too wide (13)

The measuring prism was not dried after cleaning, or the eyelids came into contact with the measuring prism whilst measuring.

➔ *The slit lamp should be withdrawn and the measuring prism dried with a cotton wool swab.*

13


#### The fluorescein ring is too narrow (14)

The lacrimal fluid has dried.

➔ *Let the patient close the eyes once or twice. Then repeat the measuring procedure.*

14


#### No semi circular rings appear, only the centre line (15)

The measuring prism does not touch the cornea! Should the patient withdraw the head slightly, then irregular pulsations will occur and the measuring prism only contacts the eye intermittently. Should the patient withdraw the head further, then the fluorescein rings will disappear altogether.

➔ *Possibly employ a headband on the patient.*

15


#### Both of the too large semi circular rings appear partly (16)

Should the slit lamp be moved forward towards the patient, or should the patient move towards the slit lamp whilst the measurement is being taken, then the feeler arm will come into contact with a sprung stop piece. The applanation surface is then too large.

➔ *The image will not change when turning the measuring drum. Withdraw the slit lamp until the regular pulsations of a corresponding smaller applanation surface indicate the correct measuring position and pressure changes lead to immediate applanation surface changes.*

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
**Position too much at left / right**

**The upper semi circular ring appears partly (17)**

The measuring prism is not focused on the eye, the eye is too far on the right.

 Move the slit lamp to the right with the guide lever.

17


**The upper semi circular ring appears completely - the lower ring partly (18)**

The measuring prism is not focused on the eye, the eye is still too far on the right.

 Move the slit lamp to the right with the guide lever.

18


**The lower semi circular ring appears completely - the upper ring partly (19)**

The measuring prism is not focused on the eye, the eye is still too far on the left.

 Move the slit lamp to the left with the guide lever.

19


**The lower semi circular ring appears partly (20)**

The measuring prism is not focused on the eye, the eye is too far on the left.

 Move the slit lamp to the left with the guide lever.


20

**Correct position! (21)**

The edges of both semi circular rings meet exactly in the centre of the ocular.

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
**Position too high / low**

**A semi circular ring appears partly in the upper half (22)**

The measuring prism is not focused on the eye, the eye is too far up.

 Move the slit lamp upwards with the guide lever.

22


**The ring appears completely in the upper half (23)**

The measuring prism is not focused on the eye, the eye is still too far up.

 Move the slit lamp upwards with the guide lever.

23


**The ring appears almost completely in the upper half, and partly cut in the lower half (24)**

The measuring prism is not focused on the eye, the eye is still too far up.

 Move the slit lamp upwards with the guide lever.

24


**Two partly cut rings appear, the larger one in the upper half (25)**

The measuring prism is not focused on the eye, the eye is still too far up.

 Move the slit lamp upwards with the guide lever.


25

**Correct position! (26)**

The edges of both semi circular rings meet exactly in the centre of the ocular

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
## Inaccurate pressure

**The outer edges of the fluorescein rings touch each other (27)**

Pressure is too strongly reduced.

 Increase the pressure by turning the measuring drum.

27


**The fluorescein rings coincide and form a line (28)**

The pressure is reduced.

 Increase the pressure by turning the measuring drum.

28


**The fluorescein rings do not touch each other (29)**

The pressure is too strongly increased.

 Reduce the pressure by turning the measuring drum.


29

**Correct position! (30)**

The inner edges of the fluorescein rings touch each other.

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### 3.6 Tips for the measurement

#### NOTICE

Only the shortest possible measuring time periods should be carried out on each eye.

Should drying of the cornea epithelium be observed, then vision and field of sight should be checked in advance.

The measuring procedure can be repeated several times. Excited or anxious patients often display a higher ocular pressure during the first measurement.

For this reason, a reduction of tension is experienced during the first few minutes, as patients notice that tonometric examinations cause no unpleasant effects. Patients will feel absolutely nothing when well anaesthetised and with the eyes wide open.

A test measuring procedure should therefore be first made on each eye. These results can be disregarded. Thereafter, three measuring procedures are to be carried out on each eye. The readings will be correct when the pressure has been stabilised. When measuring procedures have been correctly carried out, divergences in the results will be only of the order of  $\pm 0.5$  mm Hg.

When the measuring procedure lasts too long on one eye, drying will occur on the cornea epithelium of both eyes. A ring of fluorescein-positive deposits around the contact positions of the cornea and the measuring element will occur on the eye under examination.

On the other eye, map-like fluorescein-positive dried places will occur, which will not allow a reliable measuring procedure.

Extensive drying phenomena disappear after a short time without treatment becoming necessary. Focused vision will be affected by such fine epithelium defects.

### 3.7 Astigmatism

If the cornea is spherical, measurements can be made on any meridian, but it is most convenient to do it on the  $0^\circ$  meridian.

This is not so when eyes with higher corneal astigmatism than 3 dioptres are examined, as the flattened areas are not circular but elliptic.

It has been calculated that, in cases of larger corneal astigmatisms, a surface of  $7.354 \text{ mm}^2$  ( $\varnothing 3.06 \text{ mm}$ ) is to be applanated, when the measuring prism is at an angle of  $43^\circ$  to the meridian of the greatest radius.

#### Example

If the corneal astigmatism amounts to

$$6.5 \text{ mm} / 30^\circ = 52.0 \text{ D} / 30^\circ \text{ and}$$

$$8.5 \text{ mm} / 120^\circ = 40.0 \text{ D} / 120^\circ$$

the graduation value  $120^\circ$  of the prism is set at the red  $43^\circ$  mark of the prism holder.

If there is a corneal astigmatism of

$$8.5 \text{ mm} / 30^\circ = 40.0 \text{ D} / 30^\circ \text{ and}$$

$$6.5 \text{ mm} / 120^\circ = 52.0 \text{ D} / 120^\circ$$

the graduation value  $30^\circ$  is set at the red  $43^\circ$  mark. In other words, set the axial position of the greatest radius, that is the axis of a minus cylinder, on the prism graduation at the red mark on the prism holder.

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
## Appendix A Installation

**NOTICE** When mounting the tonometer onto the equipment of other manufacturers, ensure that the inter-related fitting dimensions of the various tonometer models are properly observed!

### A.1 Tonometer AT 900 C/M


A mounting base with insertion peg is required when using the Applanation Tonometer AT 900 C/M. It is to be attached as follows:

1. Remove the retaining screw on the top of the middle cylinder housing of the microscope.
2. Leave the spring in the aperture.
3. Mount the base for the tonometer and secure firmly with the screw.
4. Mount the tonometer bearer arm onto the peg of the base and swing the tonometer to the right until it engages. The tonometer also stays in this position when not in use.
5. The stop peg must be **mounted** (31) when using the BM 900® and BM 900®V, it must be **removed** (32) when using the BC 900® and BD 900®.

### A.2 Tonometer AT 900 M/Q


The Applanation Tonometer is mounted by inserting the peg of the tonometer base into one of the apertures on the horizontal guide plate over the slit lamp axis.

### A.3 Tonometer AT 900 BQ

For installation, a small screwdriver and a 5 mm Allen key is required.

1. Move the slit lamp towards the examiner up to the stop piece and block the instrument base and the microscope arm in the middle position.
2. Swing the illumination unit right round to the left.
3. Remove the breath shield from the microscope arm.
4. Remove the covering screw on the milled surface on the right of the microscope arm.

**CAUTION** The screw on the left side of the arm can be slackened with the screwdriver by unscrewing clockwise (valid only for a subsequent supply of a tonometer)!

5. Move the arm of the tonometer and the radial movement bearing together in a stretched position.
6. Support the tonometer on the guide plate and tighten the fixing screw in the boring on the microscope arm. If no guide plate is mounted, keep tonometer in a horizontal position.
7. Release the blockage on the microscope arm and the instrument base. Then install the breath shield.

### A.4 Tonometer AT 870

**NOTICE**

The instruction manuals for the equipment of other manufacturers should be observed.

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## Appendix B Maintenance

### B.1 Check tonometers AT 900

This check should be made 1 x per month. In case the check shows errors:


1. Is the measuring prism properly inserted?
2. Is the check weight precisely adjusted?
3. Repeat check procedure!

Faulty equipment is to be returned immediately to the Haag-Streit distributor.


#### NOTICE

Haag-Streit provides a repair and service package. For further information please contact us directly.


AT 900 C/M

AT 900 M/Q

AT 900 BQ


AT 870

#### Check at measuring drum setting 0


Check position - 0.05

Turn the zero calibration on the measuring drum downwards by the width of one calibration marking, against the index marker. When the feeler arm is in the free movement zone, it should then move itself against the stop piece in the direction of the examiner.

Check position + 0.05

Turn the zero calibration on the measuring drum upwards by the width of one calibration marking, against the index marker. When the feeler arm is in the free movement zone, it should then move itself against the stop piece in the direction of the patient.

#### Check at measuring drum setting 2


**CAUTION** This is the most important check procedure, as the measuring of the ocular pressure in this area is highly significant.

For this control procedure the check weight is used. Five circles are engraved on the weight bar. The middle one corresponds to drum position 0, the two immediately to the left and right to position 2 and the outer ones to position 6.

One of the marks on the weight corresponding to drum position 2 is set precisely on the index mark of the weight holder.

Holder and weight are then fitted over the axis of the tonometer (33) so that the longer part of the weight points towards the examiner.

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

Check position 1.95

When the measuring drum setting is 1.95, the feeler arm should move from the free movement zone against the stop piece towards the examiner.

Check position 2.05

When the measuring drum setting is 2.05, the feeler arm should move from the free movement zone against the stop piece in the direction of the patient.


### Check at measuring drum setting 6

Turn the weight bar to scale calibration 6, the longer part shows in the direction of the examiner.

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
Check position 6.1

The check point is 6.1. The calibration marking 6 on the measuring drum is to be turned 1/2 interval upwards. The feeler arm should move in the direction of the patient.

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


AT 870

### Check at measuring drum setting 0

Swing the feeler arm into the working position and set the measuring drum to zero.

The feeler arm with the measuring prism inserted must swing to and fro between both stop pieces, when set in motion.

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
The check weight is to be used for this control procedure.

Five circles are engraved on the weight bar. The middle one corresponds to drum position 0, the two immediately to the left and right to position 2 and the outer ones to position 6.

One of the marks on the weight corresponding to drum position 2 is set precisely on the index mark of the weight holder. Holder and weight are then fitted over the axis of the tonometer so that the longer part of the weight points towards the patient.

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

Check position 1.95

Displace the calibration marking 2 on the measuring drum downwards against the index marker by the width of the marking. When the feeler arm is slightly moved in the direction of the free movement zone, it should move independently against the stop piece in the direction of the examiner.

Check position 2.05

Displace the calibration marking 2 on the measuring drum upwards against the index marker by the width of the marking. The feeler arm should move against the stop piece on the side of the patient.


### Check at measuring drum setting 6

Turn the weight bar to scale calibration 6, the longer part shows in the direction of the patient.


Check position 5.9

The check point is 5.9. The calibration marking 6 on the measuring drum is to be turned 1/2 an interval downwards. The feeler arm should move towards the examiner.

Check position 6.1

The check point is 6.1. The calibration marking 6 on the measuring drum is to be turned 1/2 an interval upwards. The feeler arm should move in the direction of the patient.

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### B.3 Cleaning / disinfecting of measuring prisms

- Tonometer measuring prisms are not dispatched in a disinfected state and should always be disinfected before use.

- Thorough cleaning under running cold water during 30 to 60 seconds is mandatory before undertaking disinfecting routines. For more intensive cleaning in case of strong contamination carefully rub the prism with a cotton-wool ball under running water. In addition only mild soap may be used.

- For disinfecting purposes we recommend an aqueous 'Pantasept®' solution of 3% during 10 minutes. Refer to the latest versions of the manufacturer's instructions and sales literature for details of the exact concentration, effecting time and application of disinfectant.

Alternative agents: only those products which are listed on Haag-Streit AG internet homepage ([www.haag-streit.com](http://www.haag-streit.com)) may be used.

- After disinfection the measuring prisms must be rinsed thoroughly at least 10 to 60 minutes under running cold drinking water. Dry them with a clean soft one-way tissue and store them in a clean dry container.

- Only measuring prisms, free of any residue, should be used. Residue may cause provoking and caustic reactions on the patient's eye.

Not allowed!


Measuring prisms can get tension rips because of:

- Alcohol
- Acetone
- UV radiation
- Sterilization
- Immersing in fluid for more than 1 hour
- Temperature above 60 °C

### B.4 Accessories disinfecting / storing

For efficient disinfecting routines we recommend our DESINSET (9) which comprises:

- cleaning- and
- storage containers
- 'Pantasept®' disinfecting solution

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**Poids net** 0.850 kg (sans accessoires)

## Appendix C Technical specifications


### Applanation Tonometer 900

#### Measuring force generated


by leverage weight

- Installation**
- 900 C/M Mountable on peg on microscope
  - 900 M/Q Insertable in the guide plate on slewing peg for microscope and illumination unit arm
  - 900 BQ Mountable on microscope arm


AT 900 C/M

AT 900 M/Q

AT 900 BQ

AT 870

**Measuring range** 0 - 80 mm Hg

**Measurement divergence**, of the force impingement on the measuring prism over a measuring range from 0 - 58.84 mN

Standard divergence:  
 $0.49 \text{ mN} \leq 3 s \leq 1.5\% \text{ of rated value}$

**Reverse span**  $\leq 0.49 \text{ mN}$

**Net weight**

• 900 C/M	0.500 kg
(without accessories) • 900 M/Q	0.410 kg
• 900 BQ	0.680 kg

### Applanation Tonometer 870

#### Measuring force generated

by spring force

**Installation** on / over microscope

**Measuring range** 0 - 80 mm Hg

**Measurement divergence**, of the force impingement on the measuring prism over a measuring range from 0 - 58.84 mN

Standard divergence:  
 $0.49 \text{ mN} \leq 3 s \leq 1.5\% \text{ of rated value}$

**Reverse span**  $\leq 0.25 \text{ mN}$


**Net weight** 0.850 kg (without accessories)

AT 900 C/M

AT 900 M/Q

AT 900 BQ

AT 870

Certified Quality Management System

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
## **USA only:**


*Avoid errors!*

*Make the page on the right  
available to all staff mem-  
bers involved in cleaning  
and disinfecting the tonome-  
ter prisms.*

*Thank you.*

• sensitive feeler arm!

Allowed  
• mild soap

The use of the following products / procedures for disinfection will degrade the tip and cause it to fail:

- Alcohol → See page 43 for more information regarding the use of alcohol for disinfection.
- Acetone
- UV radiation
- Sterilization
- Immersing in fluid for more than 1 hour
- Temperature above 60 °C

**Maintain the tonometer**

- Clean tonometer only with a dampened cloth
- Use soap only

- Do not immerse tonometer in solution!
- No corrosive agents
- No direct spraying


**Haag-Streit USA INC**  
3535 Kings Mills Road  
Mason / Ohio 45040

Phone ++ 1 513 336 7255  
Fax ++ 1 513 336 7260  
Toll Free 800 787 5426  
eMail info@haag-streit-usa.com  
internet www.haag-streit-usa.com

See also: Appendix B.3

Recommendation: Use the Haag-Streit 'DESINSET' 2.4 Accessories (page 13)


**1**

Remove carefully measuring prism from holder

- sensitive feeler arm!


**2**

Clean **30 - 60 seconds** running, cold water

- Allowed
- mild soap

**3**

Disinfect


_____	Minutes
_____	
_____	%

The use of the following products / procedures for disinfection will degrade the tip and cause it to fail:

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
**4**


Rinse **min. 10 - max. 60 minutes** running, cold drinking water

**5**

Dry With a one-way tissue clean, soft

**6**

Store Place into container clean, dry

**Maintain the tonometer**

- Clean tonometer only with a dampened cloth
- Use soap only

- Do not immerse tonometer in solution!
- No corrosive agents
- No direct spraying

**Haag-Streit AG**

Gartenstadtstrasse 10  
 CH-3098 Koeniz, Switzerland  
 Phone ++ 41 31 978 01 11  
 Fax ++ 41 31 978 02 82  
 eMail info@haag-streit.ch  
 internet www.haag-streit.com

1500.7006000.04101

03.08 - 2