Description of unit and
operating instructions for patients

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



Contents

1. Overview	4	6. Hygienic preparation	25
2. Description of unit	8	6.1 Intervals	25
2.1 Intended use	8	6.2 Cleaning	26
2.2 User qualification	8	6.3 Disinfection, sterilization	28
2.3 Functional description	8	6.4 Change in patients	29
3. Safety instructions	10	7. Function check	30
3.1 Safety instructions	10	7.1 Intervals	30
3.2 Contraindications	12	7.2 Method	30
3.3 Side effects	13	8. Troubleshooting	34
3.4 Special labelling on unit	13	8.1 Malfunctions	34
4. Set up unit	15	8.2 Alarms	35
4.1 Set up and connect unit	15	9. Maintenance	39
4.2 Put on the Mask	16	9.1 Intervals	39
4.3 Accessories	16	9.2 Filter change	40
5. Operation	19	9.3 Change pressure measuring hose	42
5.1 Start up the VENTImotion	19	9.4 Disposal	42
5.2 Auto switch-on	20	10. Scope of supply	43
5.3 Set soft start	21	10.1 Standard scope of supply	43
5.4 Set humidifier	22	10.2 Accessories	43
5.5 After usage	23	10.3 Spare parts	45
5.6 Travelling with the VENTImotion	23	11. Technical specifications	46
		11.1 Pneumatic diagram	48
		11.2 Safety distances	49
		12. Warranty	50
		13. Index	52

1. Overview

VENTImotion

1 Bacteria filter

For protecting patient from bacteria, especially if the device is being used by several patients.

2 Connecting cord

For connecting the VENTImotion to the supply.

3 Handle

For transporting the unit.

4 Serial interfaces

Used to connect to devices for display and analysis.

5 Control panel and display

For controlling and monitoring the VENTImotion and connected accessories.

6 Connection for humidifier

For connecting the VENTIclick humidifier available as an accessory.

7 Unit outlet

Respiratory air flows out from here to the patient via the hose system and mask.

8 Sealing plug (2x)

For sealing the pressure measuring hose during cleaning.

9 Drying adapter

Required if the hose system is being dried with the VENTImotion.

10 Adapter

For connecting the hose system to the unit outlet.

11 Pressure measuring hose

For measuring the pressure in the mask and thus for controlling flow.

12 Hose system

The air flows to the mask through the hose system. The hose system consists of corrugated hose, pressure measurement tube and adapter.

13 Exhalation system

Carbon dioxide-enriched expired air escapes here during therapy.

14 Mask

Respiratory air at the necessary therapy pressure is administered to the patient via the mask.

15 Headgear

For correct and secure positioning of the mask.

16 Filter compartment cover, air inlet

For covering and securely positioning the coarse and fine dust filter.

17 Supply connection

This is where the connecting cord is attached to the unit.

18 Cable securing clip

Prevents the unit being disconnected inadvertently from the power supply.

19 Battery connection VENTIpower

For connecting the VENTIpower mobile power supply available as an accessory.

20 O₂ connection VENTI-O₂

For connecting the VENTI-O₂ oxygen connection valve available as an accessory.

21 Alarm acknowledgement switch with LED

The alarm acknowledgement switch can be used for temporary muting of alarms. The LED gives a visual indication of alarms.

22 On/off switch

For switching the VENTImotion on and off.

23 Dial

Central control for the VENTImotion, for navigating in the menu.

24 Operating switches

For rapid setting by a doctor, locked out in patient mode.

25 Menu switch

For switching between the standard display and the menu.

26 Humidifier switch with LED

For setting humidifier stage. Six stages are available. The LED indicates whether the humidifier is activated.

27 Soft start switch

Used for activating Softstart and for setting the Softstart time up to the maximum set by the physician.

28 VENTiPower

Available as an accessory, to provide mains-independent power to the VENTImotion.

29 Rating plate

Provides information about the unit, like serial number and year of manufacture.

30 VENTiClick

Available as an accessory, for humidifying and heating respiratory air.


31 Carrying bag

For transporting the VENTImotion.

32 VENTI-O₂

Obtainable as an accessory, for introducing oxygen into the mask.

Standard display during therapy:

1 Status line

This is where information about the status of the unit is displayed, like filter change or maintenance due.

2 Soft start symbol

Indicates that soft start is activated, the number showing time remaining in minutes.

3 Respiratory phase switch indicator

Indicates whether the current respiratory phase switch is spontaneous or mandatory (spontaneous: S, mandatory: T); the indicator switches from left (inspiration) to right (expiration) depending on respiratory phase; mandatory expiration is shown here.

4 Access to patient menu

The switch next to this menu item is used to switch between the patient menu and the standard display.

5 Bar chart for pressure indication

For graphical display of therapy pressure.

6 Ventilation parameters

The relevant current ventilation parameters are displayed depending on the active mode.

7 Active ventilation mode

The active ventilation mode is displayed at this point in the status line.

Symbols used in the display

Symbol	Significance
Status line:	
	Soft start active, remaining time faded in
	Filter change required
	Maintenance required
	Acoustic signal for the IPAPmin and V_{Tmin} alarms silent
	Alarm for IPAPmin and V_{Tmin} alarms deactivated
Main window:	
	Low-priority alarm triggered
	Medium-priority alarm triggered

Abbreviations used in the display

Symbol	Significance
Status line:	
S	S mode active
ST	ST mode active
T	T mode active
SX	SX mode active
SXX	SXX mode active
+v	Volume compensation activated (after mode: e.g. SXX^{+v})
Main window:	
IPAP	Inspiration pressure
EPAP	Expiration pressure
hPa	Pressure indication in Hectopascal: 1,01973 hPa equal 1 cm H ₂ O
f	Respiratory frequency
S	Spontaneously-triggered respiratory phase switch
T	Mandatorily-triggered respiratory phase switch

2. Description of unit

2.1 Intended use

The VENTImotion is a home ventilation unit for non-invasive, non-life support ventilation of adult patients with respiratory insufficiency in whom there is evidence of autonomous respiratory drive. This corresponds to the following clinical pictures:

- restrictive and obstructive ventilation disorders like paresis of the diaphragm, OSAS, COPD
- disorders of the respiratory mechanism like scoliosis, deformity of the thorax
- neurological, muscular and neuromuscular disorders
- central respiratory regulation disorders

The VENTImotion is **not suitable for life-support use**.

Use the unit only for the purpose described here.

2.2 User qualification

The VENTImotion may only be used by a doctor on an individual basis for a specific patient.

When the VENTImotion is handed over to the patient, the doctor in charge or hospital staff must instruct the patient in how the unit works.

2.3 Functional description

Provision of therapy pressure

An electronically-controlled fan draws in ambient air through a filter and delivers it to the unit outlet. From here, the air flows through the hose system and the mask to the patient.

Sensors detect the pressure in the mask and in the hose system and also the respiratory phase switch (trigger point). The fan accordingly provides the IPAP and EPAP pressures set by the doctor.

Therapy modes

The VENTImotion can be operated in six therapy modes: CPAP, S, ST, T, SX and SXX mode. Your doctor will set the mode required for your therapy on the unit.

Your doctor can activate volume compensation. Minimum volume and maximum pressure increase are set for this purpose. If the minimum volume is undershot, the unit automatically increases pressure continuously until the set maximum pressure (therapy pressure + max. pressure increase) is reached.

In time-controlled mode T and in assisted-controlled mode ST, your doctor can set respiratory frequency in the range from 6 to 45 breaths a minute and inspiration time in the range from 20 % to 67 % of the respiratory period.

In assisted modes S, SX and SXX and in assisted-controlled mode ST, your doctor can select one of 6 trigger stages for both inspiration and expiration. Your doctor can switch off the trigger for expiration. Expiration is then time-controlled.

If no breath is taken into the unit in S mode, the pressure is switched at a minimum frequency of 6 breaths a minute.

Other functions

The soft start function makes it easier to fall asleep. Your doctor sets initial pressures for inspiration and expiration which rise continuously to the therapy pressures during the soft start phase. This function can be locked out by the doctor.

The VENTImotion has an auto switch-on system. If this is activated, the unit can be switched on by a breath being taken into the mask. The unit is still switched off using the on / off switch .

The display shows therapy mode and, depending on mode, the values currently being applied for CPAP/IPAP and EPAP and for respiratory frequency (f). Spontaneous or mechanical respiratory phase switches are also displayed and the pressure change shown in the form of a graph.

3. Safety instructions

3.1 Safety instructions

Read these operating instructions carefully. They form part of the unit and must be available at all times. Use the unit only for the purpose described here (see "2.1 Intended use" on page 8).

For your own safety and for the safety of your patients and to comply with the requirements of the Medizinproduktegesetz (German law governing medical devices), please observe the following.

Operating the unit

Caution!

- Check that the supply voltage of the unit matches your supply voltage. The unit can operate with voltages of 115 V and 230 V. It sets itself to one of these voltages automatically.
- Always secure the supply plug with the cable securing clip to prevent the plug being pulled out inadvertently.
- The VENTImotion must be connected to an easily accessible supply socket so that in the event of a malfunction, the supply plug can be removed quickly.
- Do not place the unit near a source of heating or expose to direct sunlight to prevent the unit overheating. Condensation could also form in the VENTiClick humidifier which would then condense in the hosesystem.
- Do not cover the unit with blankets etc.. The air inlet would be blocked and the unit could overheat. This can lead to inadequate therapy and damage to the unit.
- Maintain a safe distance between the VENTImotion and devices which emit high-frequency radiation (e.g. mobile phones), otherwise there may be malfunctions (see page 49.)
- An alternative ventilation system should be on hand in case of a unit failure.

Caution!

- Third-party masks may be used only after approval from the manufacturer (Weinmann). Use of non-approved masks puts successful therapy in jeopardy.
- If a pneumotachograph with a high flow resistance is used in determining flow rate to start or check therapy, the trigger function of the VENTImotion may be restricted. If you have any questions, please contact the manufacturer (Weinmann).
- No antistatic or electrically-conductive hoses may be used.
- Using the VENTlick humidifier, bacteria filter and O₂ connection valve VENT-O₂ accessories may modify the characteristics of the unit. Subsequent additions to these accessories may render it necessary to readjust the unit parameters. Consult your doctor if necessary.
- Please comply with Section "6. Hygienic preparation" on page 25 to prevent infection or bacterial contamination.

Transport/accessories/spare parts/repair**Caution!**

- Do not transport the VENTImotion with the VENTlick humidifier attached. Tilting might cause water remaining in the VENTlick to run into the VENTImotion and damage it.
- Malfunctions and a lack of biocompatibility may result if third-party articles are used. Please bear in mind that in these cases, any warranty entitlement and liability shall lapse where neither the accessories recommended in the instructions for use nor original spare parts have been used.
- Have all inspections and repair work carried out only by the manufacturer (Weinmann) or specialist staff.

Oxygen supply

Warning!

- When oxygen is introduced into the respiratory flow, smoking and naked flames are prohibited. **Risk of fire.** The oxygen may be deposited in clothing, bedlinen or hair. It can only be removed by thorough ventilation.

Caution!

- Oxygen may only be introduced into the respiratory flow using the O₂ supply valve VENTI-O₂, WM 24200.
- Always follow the safety information included in the instructions for your oxygen system.

3.2 Contraindications

The VENTI*motion* should be used only with particular caution, if at all, with the following diseases. In each case, the decision whether to treat with the VENTI*motion* lies with the doctor in charge.

- Cardiac decompensation
- Severe cardiac arrhythmia
- Severe hypotension, particularly in conjunction with intravascular volume depletion
- Severe epistaxis
- High risk of barotrauma
- Pneumothorax or pneumomediastinum
- Pneumoencephalus
- Cranial trauma
- Status following brain surgery and surgical procedure performed on the pituitary gland or middle/inner ear
- Acute sinusitis, inflammation of the middle ear (otitis media) or perforated eardrum
- Dehydration


Threatening situations have not been observed with the VENTI*motion*.

3.3 Side effects


The following undesired side effects may occur during short and long-term usage of the VENTImotion:

- pressure points on the face caused by the mask and forehead padding
- reddening of skin on face
- blocked nose
- dry nose
- dry mouth in the morning
- feeling of pressure in sinuses
- irritation of conjunctiva of eyes
- gastrointestinal insufflation of air ("stomach bloating")
- nosebleeds.

3.4 Special labelling on unit

Front	
7 	Unit outlet: outlet for ambient air at 4 - 35 hPa
6 	Socket: electrical connection for the VENTi <i>click</i> humidifier, WVM 24365; max. power consumption at 40 V: 600 mA
Rear	
16 	Unit inlet: inlet for ambient air at room temperature

Side connection sockets		
4a		Connection for specialist staff to adjust therapy parameters with VENTIsupport, WM 93340; max. power consumption at 12 V: 15 mA
4b		Connection for optional additional SaO ₂ devices; max. power consumption at 12 V: 25 mA
19		Connection for supply-independent operation using VENTIpower, WM 24100; power supply: 12 V/40 V DC
20		Connection for controlling the O ₂ connection valve VENTI-O ₂ , WM 24200. Max. power consumption at 12 V: 125 mA
Rating plate		
	Do not dispose of the unit in domestic waste!	
	Protection class B	
	Protection class II, protective insulation	
	Year of manufacture	


4. Set up unit

4.1 Set up and connect unit

Set up the unit on a flat surface, e.g. on a bedside cabinet or on the floor next to the bed. Keep a space of min. 5 cm between the rear of the unit and the wall, as the air inlet is on the rear of the unit.

Caution!

Do not cover the unit with blankets etc.. The air inlet would be blocked and the unit could overheat. This can lead to inadequate therapy and damage to the unit.


1. Connect the connecting cord to the supply connection of the unit. To do so, lift the cable securing clip, plug the plug into the supply connection and clip the cable securing clip over the plug.

Caution!

Always secure the supply plug with the cable securing clip to prevent the plug being pulled out inadvertently.

2. Connect the connecting cord to a supply socket. The VENTImotion is designed for a supply voltage of 115 V ~ to 230 V ~. "Weinmann VENTImotion" appears in the display.
3. Plug the hose system adapter in the unit outlet. The unit is now operational.

4.2 Put on the Mask

The VENTImotion is designed for use with nasal, nasal/oral and full facial masks. Proceed as follows to put on the mask.

1. Adjust the forehead support of the I mask (if provided).
2. Connect the headgear/the headband to the mask and position the mask.
3. Adjust the headgear/headband so that the bulge of the mask only results in low pressure to prevent pressure points on the face.

Please refer to the operating instructions for the relevant mask.


4.3 Accessories

Separate exhalation system

A separate exhalation system is required if you are using a mask without an integrated exhalation system.

Used carbon dioxide (CO₂)-enriched air escapes from the mask through the exhalation system. Without an exhalation system, the CO₂ concentration would rise to critical values in the mask and hose, thus hindering your breathing.


The exhalation system also allows you to breathe through the nose or mouth if the unit malfunctions. In the event of a malfunction with a full facial mask, you breathe through a safety valve on the mask.

The sleeve on the end of the hose system has a diameter of 19.5 mm and fits over a 22 mm standard cone. Put the exhalation system on the sleeve of the creased hose.


Please refer to the operating instructions for your exhalation system.

Humidifier VENTi**click**

The VENTi**click** humidifier is positioned between the VENTI**motion** and the hose system. The inlet connection stub and the electrical connection of the heating element must be pointing at the VENTI**motion**. Please also see the operating instructions for the VENTi**click**.

Oxygen supply using the VENTI-**O₂**

The supply of oxygen must have been prescribed by the doctor in charge.


For safety reasons (risk of fire), the direct introduction of oxygen into the hose system or the mask is not permitted without a special protective device.

On the VENTI**motion**, the introduction of oxygen is permitted solely using the VENTI-**O₂** oxygen valve, WM 24200. This valve allows up to 4 l/min. of oxygen to be introduced. In the event of a fault, the VENTI-**O₂** gives off the oxygen into the ambient air. This means that it cannot accumulate in the unit.

The oxygen can be supplied using an oxygen concentrator (e.g. Oxymat 3), central gas supply system, liquid oxygen or an oxygen cylinder with a suitable pressure reducer. The external oxygen source must be equipped with a flow-adjustment device independent of the VENTI-**O₂**.

You must follow the safety instructions for handling oxygen without fail, as well as the operating instructions for the VENTI-**O₂** oxygen valve and the oxygen device used.

Bacteria filter


If the *VENTImotion* is intended for use by several patients (e.g. in hospital), bacteria filter WM 24148 should be used to protect against infection. It is installed between the hose system and the *VENTImotion* or the *VENTIclick*. Please also observe the enclosed operating instructions for the bacteria filter.

Caution!


The bacteria filter represents an additional resistance in the air stream. This may lead to a change in the response characteristics of the trigger. For this reason, you should have the unit parameters reset if a bacteria filter is fitted retrospectively.

VENTIpower

For mains-power independent operation *VENTIpower* is available.

VENTIpower can also be connected to *VENTImotion* in parallel to mains power. *VENTIpower* then remains in standby mode. If mains power fails, *VENTIpower* automatically supplies power to the therapy device with a delay of about 4 seconds.


Using *VENTIpower* you can operate *VENTImotion* independent from mains power in the carrying bag WM 24888. To do so place *VENTIpower* in the dedicated section of the carrying bag and connect it to the therapy device.

Please also observe the operating instructions for *VENTIpower*.

5. Operation

5.1 Start up the VENTImotion

1. If your mask has no integrated exhalation system, place the exhalation system on the end of the hose system (see "Separate exhalation system" on page 16).

Caution!

Always use an exhalation system. Otherwise the CO₂ concentration would rise to critical values in the mask and hose, thus hindering your breathing.


2. Put on the mask (see "4.2 Put on the Mask" on page 16) and connect the hose system including the exhalation system to the mask (plug cone). Please follow the relevant operating instructions for the mask and exhalation systems.


3. Press the on / off switch . If auto switch-on is activated, you can also position the mask and switch on the VENTImotion by taking a breath (see "5.2 Auto switch-on" on page 20).

Operating hours and the Weinmann software version will appear in the display for approx. 3 seconds. The unit starts to convey air through the hose system. The display switches to the standard display.

The display shows the therapy mode set, therapy pressures (only CPAP pressure in CPAP mode) and current respiratory frequency.

When soft start is activated, the soft start indicator  appears with the remaining soft start time, max. 30 minutes, or the max soft start time specified by the doctor (see "5.3 Set soft start" on page 21).

The bar chart shows the inspiratory and expiratory pressure curve.


The respiratory phase switch indicator shows whether the current respiratory phase was triggered spontaneously by the patient (**S**) or mechanically (**T**). Depending on the respiratory phase, the indicator switches from left (inspiration) to right (expiration). The illustration shows mechanically-triggered expiration.

5.2 Auto switch-on


When auto switch-on is activated, you can switch on the VENTImotion using the on / off switch or by "taking a breath". You can only switch off the VENTImotion using the on / off switch .

Activate auto switch-on


1. Seal the hose system or put on the mask.
2. Press the on / off switch to switch on the VENTImotion.


3. Press the menu switch . The patient menu appears in the display. The selection bar is over **Drying process**.

4. Use the dial to select the menu item **Auto switch-on (Off)** and press the dial.

The message **Auto switch-on ON** appears for three seconds. After the message, the VENTImotion reverts to the patient menu. Auto switch-on is now activated.

- Exit the patient menu using the dial (move selection bar to **Back** and press dial) or using the menu switch .


If no switch is activated for 5 minutes, the display reverts to the standard display.

Deactivate auto switch-on

Repeat steps 2.–5. as described in „Activate auto switch-on“. When the automatic system is switched off, you can only switch on the VENTImotion using the  switch.

5.3 Set soft start

If your doctor has enabled the soft start function, the soft start time can be selected in 5-minute increments up to a maximum time of 30 minutes. Your doctor can limit the maximum time to less than 30 minutes. Proceed as follows to set soft start time.

- Seal the hose system or put on the mask.
- Switch on the VENTImotion and hold down the soft start switch  until the soft start window appears. Change the soft start time using the dial (turn to the right to increase time and turn to the left to decrease time). Alternatively, keep pressing the soft start switch  to increase the soft start time in 5-minute increments (once maximum time is reached, the unit restarts at the minimum time of 5 minutes if pressed again).


- Press the menu switch  or the dial to store the soft start time. The soft start time displayed will be stored and the soft start window closes automatically. If no switch is pressed for 4 seconds, the soft start time displayed is likewise stored. The settings are retained after the unit is switched off.

Note

VENTI*motion* starts in soft start mode automatically if this was activated when the unit was last used. Soft start can be switched on and off at any time by briefly pressing the soft start switch .

5.4 Set humidifier

The VENT*click* humidifier can be used to humidify and heat up the respiratory air supplied by the VENT*motion*. Heat output can be selected in six stages. Proceed as follows.

1. Seal the hose system or put on the mask.
2. Switch on the VENT*motion* and press the humidifier switch  until the Humidifier stage window appears. Alter the heat output (stages 1-6) of the humidifier using the dial or repeatedly press the  switch.
3. Press the menu switch  or the dial to store heating output. The Humidifier stage window closes automatically. If no switch is pressed for 4 seconds, the heating output displayed is likewise stored. The settings are retained after the unit is switched off.


Please also see the operating instructions for the VENT*click*.

Note

- VENT*motion* automatically starts with the VENT*click* activated if this was activated when the unit was last used. VENT*click* can be switched on and off at any time by briefly pressing the Humidifier switch . When the humidifier is activated, the green status indicator next to the humidifier switch is illuminated.
- If you press the humidifier switch  without a humidifier being connected, the VENT*motion* switches this function off again automatically.

5.5 After usage

1. Remove the headgear/headband and the mask.
2. Keep the on / off switch  depressed for 2 seconds to switch off the unit. The fan switches off. Date and duration of the last therapy session appear in the display. "Weinmann VENTImotion" appears in the display. The unit then switches to standby mode. "Standby" appears in the display.


3. Remove the hose connection and, if necessary, the exhalation system from the mask.
4. Clean the mask and the exhalation system (see "6. Hygienic preparation" on Page 25).

Note!

To save electricity, you can remove the plug of the supply connection from the supply socket. The stored therapy parameters and settings will be retained.

5.6 Travelling with the VENTImotion


VENTImotion should only be transported long distances in the carrying bag provided.

Place the unit, the connecting cord, the hose system including the drying adapter, the mask including the exhalation system and, if necessary, the VENTlick humidifier and the VENTI-O₂ oxygen valve in the carrying bag. Please do not forget to take spare filters and all the operating instructions with you.

Warning!

- Never operate the VENTI-O₂ oxygen valve available as an accessory in the carrying bag. **Risk of fire.** Oxygen may be deposited in the carrying bag.

If you want to use the VENTImotion on the journey, leave the unit and the VENTIpowers in the bag during operation to prevent damage. To operate the unit in the bag, see the operating instructions for the VENTIpowers.

If you want to take the VENTImotion on an aircraft as hand baggage, find out about any formalities required before departure. You can obtain a certificate for transportation in an aircraft from the manufacturer (Weinmann).

6. Hygienic preparation

6.1 Intervals


You should check the filters and wipe clean the housing and the filter compartment cover regularly. The headgear/headband can be washed as necessary. Please also observe the following intervals.

Interval	Activity
Daily	<ul style="list-style-type: none">– Clean mask– Clean hose system– Clean bacteria filter in accordance with the operating instructions– Clean VENTlick humidifier in accordance with the relevant operating instructions.– Clean the exhalation system after every use in accordance with the operating instructions .
Every 24 operating hours	<ul style="list-style-type: none">– Change the particle filter in the bacteria filter
Weekly	<ul style="list-style-type: none">– Clean the mask thoroughly in accordance with the operating instructions.– Clean coarse dust filter
Every 1000 operating hours	<ul style="list-style-type: none">– Change fine filter (filter change indicator ). If it is soiled, change it sooner
Every 6 months	<ul style="list-style-type: none">– Change coarse dust filter. Change it sooner if soiled or worn– Change pressure measuring hose – sooner if soiled – (see “9.3 Change pressure measuring hose” on page 42).
Annually	<ul style="list-style-type: none">– Change hose system.

6.2 Cleaning

Clean hose system

1. Detach the hose system from the unit and the exhalation system.

2. Remove one end of the pressure measuring hose (shake gently if necessary) and seal it with the sealing plug provided. At the other end, use the second sealing plug to seal the small opening of the adapter to prevent water getting inside.

3. Clean the creased hose so no residues remain using a small amount of detergent in warm water. Rinse the inside of the hose thoroughly in the process.


4. Rinse both the inside and the outside of the creased hose thoroughly using clean warm water.

5. Shake the hose system thoroughly.

6. Hang up the hose system and allow to drip to prevent moisture getting into the VENTImotion.

7. Remove the plugs from the pressure measuring hose.


Dry the hose system using the VENTImotion

1. If water accidentally gets into the pressure measuring hose, insert the red drying adapter supplied in the unit outlet.

2. Insert the adapter of the hose system in the unit outlet or if necessary, onto the red drying adapter.

3. To activate the drying process, press the on / off switch . Use the menu switch  to reach the patient menu. Press the dial to start the drying process. The unit now dries the hose system.

The remaining drying time is displayed. The VENTImotion switches off automatically after 30 minutes. You can interrupt the process at any time by keeping the menu switch , the dial or the on / off switch  depressed for 2 seconds. If the hose system is still damp in places after drying, start the drying process again.

4. Remove the drying adapter from the unit outlet.


Clean the VENTImotion housing

Warning!

- There is a risk of electric shock. Before cleaning, remove the connecting cord from the supply connection on the rear of the unit and the supply plug from the supply socket without fail.
- Make sure no liquids get into the unit. Never immerse the unit in disinfectant or other liquids. Otherwise the unit may be damaged, resulting in a hazard to users and patients.

1. Wipe the unit and the connecting cord clean with a soft damp cloth. The VENTImotion must be completely dry before being started up.

2. Remove the filter compartment cover.
 3. Remove the coarse dust filter as described in 9.2.
 4. Clean the filter compartment cover under running water so no residues remain. Then dry it carefully.
5. Replace the coarse dust filter and the filter compartment cover as described in 9.2.

Clean coarse dust filter/change fine filter

1. Remove the filter compartment cover as described in 9.2.
2. Remove the coarse dust filter from the filter compartment cover and clean it under clear running water so no residues remain.
3. Change the fine filter if required.

4. Let the coarse dust filter dry. The coarse dust filter must be completely dry before the unit is started up.
5. Replace the coarse dust filter and close the filter compartment cover.

Note:

the fine filter cannot be cleaned. It is changed every 1000 operating hours.

Accessories

To clean the exhalation system, the mask, the headgear/headband, the VENTIpower, the VENTI-O₂, the bacteria filter and the VENTIconclick, please see the section entitled "Cleaning" in the relevant operating instructions.

6.3 Disinfection, sterilization

If necessary, for example after an infectious illness or if the unit has become unusually soiled, you can also disinfect the housing, the connecting cord, the hose system, the bacteria filter housing, the VENTIconclick, the VENTI-O₂ and the VENTIpower. Refer to the operating instructions for the disinfectant used. We recommend that you wear suitable gloves (e.g. household or disposable gloves) when disinfecting the equipment.

Unit

The housing and the connecting cord of the VENTI-motion are cleaned by simply being wiped over with disinfectant. We recommend TERRALIN for this purpose.

Hose system

We recommend GIGASEPT FF as disinfectant. When using GIGASEPT FF, use the same procedure as described under "Cleaning".

After disinfecting, rinse all parts thoroughly in distilled water. Let the parts dry completely.

Let the hose system drip. Dry the hose system using the VENTI-motion as described in 6.2.

- Creased hose WM 241 30 (clear) can be cleaned using hot water up to 70 °C. Sterilization is not permitted.

- Creased hose WM 24120 (grey) can be steam-sterilized using equipment to EN 285. Temperature: 134 °C, minimum treatment time 3 minutes. Follow EN 554/ ISO 11134 with regard to validation and monitoring.

Accessories

To disinfect the exhalation system, the mask, the VENTIpower, the VENTI-O₂, the bacteria filter and the VENTlick, please see the section entitled "Disinfection and sterilization" in the relevant operating instructions.

6.4 Change in patients

If the unit is being operated with a bacteria filter, please note the following:

- change the bacteria filter

or

- sterilize the bacteria filter and change the particle filter located inside it.

If the unit is to be utilized for another patient without using a bacteria filter, it must undergo hygienic preparation beforehand. This has to be carried out by the manufacturer (Weinmann) or a dealer.

The procedure used for hygienic preparation is described in the inspection sheet and in the service and repair instructions for the VENTI*motion*.

7. Function check

7.1 Intervals

Perform a function check at least every 6 months. If you discover defects during the function check, you must not use the VENTImotion again until these defects have been rectified.

7.2 Method

1. Assemble the VENTImotion so that it is ready to function with the hose system, exhalation system and the supply cable connection.
2. Close the opening of the hose system, for example using your thumb or hand.
3. Switch on the unit by pressing the on / off switch .
4. If soft start is switched on, switch it off by pressing the soft start switch .

Depending on the operating mode set, now check the following functions:

Function: / Mode:	CPAP	S	T	ST	SX	SXX
Pressure setting	•	•	•	•	•	•
Triggering	–	•	–	•	•	•
Soft start	•	•	•	•	•	•
Flow sensor	–	•	•	•	•	•
Alarms	•	•	•	•	•	•

If the values/functions set out below are not met/fulfilled, send the unit to your dealer or to the manufacturer (Weinmann) for repair.

Pressure setting

1. Wait approx. 1 minute, then read the CPAP pressure displayed or the IPAP and EPAP pressures displayed in the bar chart and check whether the values displayed correspond to the prescribed values.

Triggering

Note

This check is not necessary when the unit is operated in CPAP or T mode.

1. Switch off the unit by keeping the on / off switch  depressed for 2 seconds.
2. Fit the mask.
3. Switch on the unit by pressing the on / off switch  .
4. Put on the mask and breathe normally. In ST mode, your respiratory frequency must be above the prescribed frequency at this point so that the unit's time-controlled trigger is not activated.

The unit must react to the respiratory phase switch by switching pressure level .

Note

If the trigger for expiration does not function, it may have been switched off. Ask your doctor whether this is the case.

Soft start:

Note


This test is not necessary if soft start has been locked out by the doctor.

Switch on soft start by pressing the soft start switch .


The soft start indicator  appears and the soft start time is displayed.

Flow sensor/pressure sensor


1. Plug the red drying adapter supplied into the unit outlet.


2. To start the function check for the flow sensor, press the on / off switch . Use the menu switch  to reach the patient menu.


3. Use the dial to move the black selector bar to **Function check** and press the dial.

The message window "Function check running!" opens. The remaining duration of the function check is displayed.

4. Once the function check has been completed successfully, the message "Function check OK!" appears and the VENTImotion returns to the standard display.

5. If an implausibility is detected during the function check, the message "Flow sensor implausible!", "Sensor system implausible!" or "Implausible pressure measurement" appears.

Close the window using menu switch , switch off the VENTImotion and disconnect the unit from the supply for 5 minutes. Check whether the drying adapter has been attached correctly. Perform the function check again.

If the messages "Flow sensor implausible!", "Sensor system implausible!" or "Implausible pressure measurement" appear again, contact the manufacturer (Weinmann) or your dealer immediately to have the unit repaired.

Alarms

The VENTImotion performs a self-test of the sensor system when the on / off switch  is pressed, this checks the function of the alarm system. If an error occurs in the course of the self-test, an error message appears on the main screen (see also "8. Troubleshooting" on Page 34).

1. Check the buzzer and the status displays:

Every time the VENTImotion is switched on, check that the buzzer sounds and that the status display is illuminated.

2. Check the supply failure alarm:

Start the VENTImotion. Now take the supply cable out of the supply socket, the display vanishes and the buzzer sounds.

Note: the VENTImotion must have been connected to the supply for at least 5 minutes before the test is performed.

Switch the unit off again after the function check is complete by pressing the on / off switch .

8. Troubleshooting

8.1 Malfunctions

Malfunction/Error message	Cause	Rectification
No running noise, no Display	No supply voltage available.	Check connecting cord is securely connected. Possibly check function of supply socket by connecting another piece of equipment (e.g. a lamp).
Unit cannot be switched on by taking a breath.	Auto switch-on not activated.	Activate auto switch-on (5.2, Page 20).
Soft start will not switch on	Soft start function is locked out.	Check with your doctor whether the soft start function can be enabled for your therapy.
Soft start time cannot be set to 30 min.	Maximum soft start time limited by doctor.	Check with your doctor whether the maximum soft start time for your therapy can be extended to 30 minutes.
Unit runs but does not reach the set IPAP pressure.	Filter soiled.	Clean/change filter (9.2, Page 40)
	Mask leaking.	Adjust headgear/headband so that the Mask does not leak - replace if necessary
	Mask defective.	Replace mask
 Change filter!	Filter soiled.	Press alarm acknowledgement switch Clean/change filter (9.2, Page 40)
Filter change indicator  illuminated		Clean/change filter (9.2, Page 40)
 Battery discharged	Internal battery empty.	Press alarm acknowledgement switch, Battery must be changed by Weinmann or a dealer to make sure that therapy is recorded correctly.

Malfunction/Error message	Cause	Rectification
 Clock not set	Internal clock of VENTImotion no set.	Press alarm acknowledgement switch, clock must be set by your doctor to make sure therapy is recorded correctly.
 Service indicator illuminated		Unit must be serviced by Weinmann or a dealer as soon as possible.

If there are any faults which cannot be rectified immediately, contact the manufacturer (Weinmann) or your dealer immediately to have the unit repaired. Do not continue using the unit to avoid causing greater damage.

8.2 Alarms

A distinction is made between two types of alarm on the VENTImotion:

- low-priority alarms, indicated by the symbol  in the alarm window, a **continuous** yellow status indicator and an acoustic alarm (buzzer).
- medium-priority alarms, indicated by the symbol   in the alarm window, a **flashing** yellow status indicator and an acoustic alarm (buzzer).

There are no so-called "high-priority alarms" on the VENTImotion, as this unit may not be used for life-support ventilation.

Silence alarms

The acoustic alarm of the physiological alarms $V_{T\min}$ and $IPAP_{\min}$ can be deactivated by your doctor ( symbol in the status line). In this event, only the relevant alarm message appears in the display and the yellow status indicator is continuously on.



Deactivate alarms

If the symbol  appears in the status line, the physiological alarms $V_{T\min}$ and $IPAP_{\min}$ have been deactivated by your doctor.

Note

For questions about deactivating/silencing the physiological alarms, please consult your doctor.

Acknowledge alarms

If an alarm is triggered by a fault (in this case: disconnection alarm), press the alarm acknowledgement switch . The acoustic alarm will pause for approx. 120 seconds.

Once the acoustic alarm has been acknowledged, the standard display reappears. The fault which has not yet been rectified continues to be displayed in the status line and the yellow status indicator flashes (or is illuminated) until the malfunction has been rectified.

If the fault is not rectified within 120 seconds of the acknowledgement, the acoustic alarm (buzzer) sounds again.

You will find assistance in rectifying faults in the tables which follow.

Physiological alarms

Indicator	Alarm	Cause	Rectification
 	Minimum respiratory volume set by the doctor not reached. Low priority	Filter soiled.	Clean/change filter
		Mask leaking.	Adjust headgear/headband so that the mask does not leak - replace if necessary
		Mask defective.	Replace mask
		Implausible settings.	Have settings checked by your doctor.

Indicator	Alarm	Cause	Rectification
IPAPmin (leakage) 	Minimum therapy pressure set by doctor not reached. Low priority	Filter soiled.	Clean/change filter
		Mask leaking.	Adjust headgear/headband so that the mask does not leak - replace if necessary
		Masal mask defective.	Replace mask
		Implausible settings	Have settings checked by your doctor.

Technical alarms

Indicator	Alarm	Cause	Rectification
Pressure measuring hose disconnected! 	Low priority	Pressure measuring hose not properly fitted or forgotten during assembly	Check correct fitting/presence of pressure measuring hose
Device Error!   Excessive pressure!	Medium priority	Pressure sensor defective	Have unit repaired
	Pressure measuring hose blocked. Medium priority	Drops of water in pressure measuring hose	Dry pressure measuring hose as described in 6.2
Disconnection!  	Medium priority	Hose system not connected to unit properly or not connected at all	Check connection of hose to unit
		Unit being operated with open mask (mask not on).	Put on mask or switch off unit.
Device Error!   Pressure sensor implausible	Medium priority	Drops of water in pressure measuring hose	Dry pressure measuring hose as described in 6.2

Indicator	Alarm	Cause	Rectification
Device Error!  Upper temperature limit exceeded	Medium priority	Unit overheating due to direct sunlight	Allow unit to cool down, find a more suitable place to set it up.
	Medium priority	Unit operated outside permissible temperature range	
Display gone out	Acoustic signal for at least 120 seconds, no display visible. Medium priority	No supply voltage available.	Check connecting cord is securely connected. Possibly check function of supply socket by connecting another piece of equipment (e.g. a lamp).
		The VENTIpower battery has discharged.	The VENTIpower must be disconnected from the VENTImotion and recharged. The VENTImotion can continue to be operated via a supply socket.
Device Error!  Code	Medium priority	Problems with electronics or program "crashed".	Disconnect connecting cord and reconnect. Repeat after a short period if necessary.
		Sealing plugs still in Pressure measuring hose after hygienic preparation	Remove sealing plugs. Remove and reconnect connecting cord. Switch unit on again without hose system connected.
		Drops of water in pressure measuring hose	Dry pressure measuring hose as described in 6.2

If there are any faults which cannot be rectified immediately, contact the manufacturer (Weinmann) or your dealer immediately to have the unit repaired. Do not continue using the unit to avoid causing greater damage.

9. Maintenance

9.1 Intervals

We recommend only having servicing and maintenance performed by the manufacturer (Weinmann) or authorized dealers.

Both filters must be regularly checked for soiling

- The coarse dust filter should be cleaned once a week and changed no less than every 6 months.
- The fine filter should be changed after no more than 1000 operating hours (filter change symbol  appears in the display).

For reasons of hygiene, we recommend replacing the following parts at the intervals stated:

- pressure measuring hose every 6 months and earlier if soiled
- complete mask system every 6 to 12 months depending on soiling
- exhalation system in accordance with the relevant operating instructions

In addition, servicing should be carried out at the following intervals as a preventive measure:

- after every 5000 operating hours (maintenance symbol  appears in the display)
- after no more than 2 years (see servicing sticker on rear of unit)

Servicing includes:


- filter change
- cleaning unit
- replacing any defective parts
- a complete check of unit functions and pressure indicators

9.2 Filter change


Change coarse dust filter

Use only original filters made by Weinmann. If third-party filters are used, this will invalidate the warranty and may cause restrictions in terms of function and biocompatibility.

1. If the VENTi**click** humidifier is connected, first disconnect it from the unit. This will prevent water running into the unit when the filter is changed. Please also observe the operating instructions enclosed.


2. Press on the latch of the filter compartment cover and lift it off.
3. Take the coarse dust filter out of the filter compartment cover.

4. Insert the clean coarse dust filter in the filter compartment cover.
5. Insert the filter compartment cover into the opening in the housing bottom edge first.
6. Push the filter compartment cover into the housing until the latch engages.


Change fine filter

The fine filter needs changing when it has turned dark in colour, but in any event, after no more than 1000 operating hours. In the latter case, the message "Change filter!" appears in the display if the VENTi**motion** is switched on.

Press the alarm acknowledgement button . The filter change symbol  then appears continuously in the status line.


To change the fine filter proceed as follows.


1. Press the flap of the filter compartment cover. Remove the filter compartment cover.
2. Replace the fine filter with a new fine filter WM 15026. Dispose of the old filter with domestic waste.
3. Insert the lower edge of the filter compartment cover into the opening of the housing. Then push the filter compartment cover into the housing until it clicks into place.

Reset filter change indicator

After you have changed the fine filter you must reset the filter change indicator. This is necessary even if the filter was changed before the 1000 operating hours expired so the filter change symbol is not in the display.


1. To reset the filter change indicator, press the menu switch  with the unit switched on and use the dial to select the item **Filter change** from the Patient menu. Push the dial to call up the "Filter change" menu.

2. The question "Reset filter change?" appears. Use the dial to select **Yes** and confirm the selection by pressing the dial.

If you wish to cancel the process, use the dial to select **No** and press the dial. The process is aborted.


Once you have selected and confirmed **Yes** with the dial, the message "Filter change reset!" appears for about 3 seconds.


Bacteria filter

If bacteria filter WM 24148 is used, change the particle filter in the bacteria filter in accordance with the relevant operating instructions.

9.3 Change pressure measuring hose

1. Detach the sleeve on the creased hose from the adapter .
2. Remove the pressure measuring hose from the creased hose.
3. Remove the pressure measuring hose from the adapter.
4. Fit the new pressure measuring hose to the adapter.
5. Hold up the creased hose and insert the free end of the new pressure measuring hose.
6. Slide the sleeve of the creased hose onto the adapter.
7. Dispose of the old pressure measuring hose.

9.4 Disposal

Device

Do not dispose of the unit with domestic waste. To dispose of the unit properly, please contact a licensed, certified electronic scrap disposal merchant. This address is available from your Environment Officer or from your local authority.

10. Scope of supply

10.1 Standard scope of supply

VENTImotion

WM 24800

Parts	Order number
VENTImotionBasic unit	WM 24810
Hose system	WM 24130
Drying adapter	WM 24203
Carrying bag for the VENTImotion	WM 24888
Connecting cord	WM 24177
Coarse dust filter	WM 24880
Fine filter, packed	WM 15026
Operating instructions for the VENTImotion	WM 16838
Patient record	WM 16763

10.2 Accessories

The following accessories are not included in the scope of supply and must be ordered separately if required.

Parts	Order number
Humidifier VENTlick	WM 24365
VENTIsupport evaluation software (only for medical and technical staff)	WM 93340
Hose system, sterilizable, consisting of:	WM 24120
– creased hose, sterilizable	WM 24122
– pressure measuring hose	WM 24038
– hose adapter	WM 24123
– pressure connector	WM 24129
- sealing plug (2x)	WM 24115
Bacteria filter, complete	WM 24148

Parts	Order number
VENTI-O ₂ Oxygen valve	WM 24200
VENTIpower battery incl. case	WM 24100
SOMNOmask, size S	WM 25110
SOMNOmask, size M	WM 25120
SOMNOmask, size L	WM 25130
SOMNOmask, size S, blue	WM 25210
SOMNOmask, size M, blue	WM 25220
SOMNOmask, size L, blue	WM 25230
SOMNOplus, size S	WM 25710
SOMNOplus, size M	WM 25720
SOMNOplus, size L	WM 25730
Silentflow	WM 23600
Noise insulation, complete (alternative exhalation system)	WM 23685
SOMNOcap, size M, blue	WM 23673
SOMNOcap, size L, blue	WM 23674
SOMNOcap, size M, white	WM 23686
SOMNOcap, size L, white	WM 23687
SOMNOstrap, blue	WM 23695
SOMNOstrap, white	WM 23177
SOMNOstrap, extra	WM 25250

When using other mask systems, please follow the associated operating instructions.

10.3 Spare parts

Parts	Order number
Hose system, consisting of:	WM 24130
– creased hose, can be disinfected	WM 24108
– pressure measuring hose	WM 24038
– hose adapter	WM 24123
– pressure connector	WM 24129
- sealing plug (2x)	WM 24115
Coarse dust filter	WM 24880
Fine filter, packed	WM 15026
Set of spare filters necessary for one year, packed (3x fine filter, 2x coarse filter)	WM 15682
Carrying bag for the <i>VENTImotion</i>	WM 24888
Connecting cord	WM 24177
Drying adapter	WM 24203

11. Technical specifications


	VENTImotion	VENTImotion with VENTiClick
Product class as per EC Directive 93/42/EEC	IIa	
Dimensions WxHxD in cm	23 x 12.5 x 34	23 x 12.5 x 45.5
Weight	approx. 4.5 kg	approx. 4.8 kg
Temperature range – operation – storage	+5 °C to +35 °C –40 °C to +70 °C	
Air pressure range	600 – 1100 hPa	
Electrical connection	115 – 230 V AC, 50–60 Hz Tolerance -20 % +10 %	
Power consumption during – operation – standby	230 V 115 V 0.2 A/0.4 A 0.043 A/0.068 A	230 V 115 V 0.35 A/0.7 A 0.043 A/0.068 A
Classification as per EN 60601-1 – type of protection from electric shock – level of protection from electric shock	Protection class II Type B	
Electromagnetic compatibility (EMC) as per EN 60601-1-2 – radio interference suppression – immunity to interference	EN 55011 EN 61000-3-2, EN 61000-3-3, EN 61000-4-2 to 6, EN 61000-4-8, EN 61000-4-11	
Mean sound pressure level/ operation as per EN ISO 17510 at a distance of 1 m from unit in patient position	approx. 35 dB (A) at 20-35 hPa approx. 31 dB (A) at 20 hPa approx. 29 dB (A) at 15 hPa approx. 27 dB (A) at 12 hPa approx. 25 dB (A) at 10 hPa approx. 23 dB (A) at 7 hPa	
Sound pressure level, alarm	at least 45 dB (A)	

	VENTImotion	VENTImotion with VENTiClick
IPAP pressure range EPAP pressure range CPAP pressure range Pressure accuracy Increment	6 to 35 hPa 4 to 20 hPa 4 to 20 hPa ±0.6 hPa 0.2 hPa (1 hPa = 1 mbar ≈ 1 cm H ₂ O)	
Minimum stable limit pressure (PLS _{min}) (min. pressure in event of a fault) Maximum stable limit pressure (PLS _{max}) (max. pressure in the event of a fault)	≥ 0 hPa ≤ 60 hPa	
Respiratory frequency Accuracy Increment	6 to 45 l/min ± 0.5 l/min 1 l/min	
I:E (Ti/T): Inspiration time Increment Accuracy	20 % to 67 % of respiratory period 1 % ±1 %	
Trigger stage	can be adjusted in 6 stages, separately for inspiration and expiration. Trigger for expiration can be switched off in ST mode	
Pressure increase speed	can be set in 6 stages	
Pressure drop speed	can be set in 6 stages	
Accuracy of volume measurement	at 23 °C: ± 15 %	
Flow at max. speed at: 0 hPa Tolerance	300 l/min ±15 l/min	275 l/min ±15 l/min
Flow at max. speed with bacteria filter at: 0 hPa Tolerance	270 l/min ±15 l/min	250 l/min ±15 l/min
Heating of respiratory air as per HMV	2.5 °C	depending on heating stage
Pressure constant measured to DIN EN ISO 17510 in CPAP mode	at 20 hPa: Δp ≤ 1 hPa at 14 hPa: Δp ≤ 1 hPa at 10 hPa: Δp ≤ 1 hPa at 7 hPa: Δp ≤ 0.5 hPa	
Fine filter separation level to 2 μm	≥ 99.7 %	

	VENTImotion	VENTImotion with VENTiClick
Fine filter service life	1000 hours with normal room air	
Permitted humidity for operation and storage	≤ 95 % rH (no condensate formed)	
System resistance at an air flow of 60 l/min at the patient connection opening	VENTImotion with hose system WM 24130 and Silentflow WM 23600	VENTImotion with O ₂ hose system WM 23737, VENTiClick WM 24365 and bacteria filter WM 24148
	0.19 $\frac{\text{kPa} \cdot \text{s}}{\text{l}}$	0.29 $\frac{\text{kPa} \cdot \text{s}}{\text{l}}$

11.1 Pneumatic diagram

11.2 Safety distances

Recommended safe distances between portable and mobile high-frequency telecommunication devices (e.g. mobile phones) and the VENTImotion			
Rated output of HF device in W	Safe distance as a function of transmission frequency in m		
	150 KHz - 80 MHz	80 MHz – 800 MHz	800 MHz – 2.5 GHz
0.01	0.12	0.12	0.24
0.1	0.37	0.37	0.74
1	1.17	1.17	2.34
10	3.69	3.69	7.38
100	11.67	11.67	23.33

Contact Weinmann for further technical specifications. More details can also be found in the hospital manual and in the service and repair guide.

The right to make design modifications is reserved.

12. Warranty

- Weinmann offers a warranty that the product, when used in accordance with requirements, will remain free from defects for a period of two years from date of purchase. For products whose durability is clearly indicated as less than two years, the warranty expires on the expiration date indicated on the packaging or in the user's manual.
- Claims against the warranty can be made only when accompanied by the sales receipt, which must show salesperson and date of purchase.
- We offer no warranty in the case of:
 - Disregard of usage instructions
 - Operating errors
 - Improper use or handling
 - failure to use original spare parts
 - Third-party intervention by unauthorized persons for the purpose of device repair
 - Acts of God, e.g., lightning strikes, etc.
 - Transport damage as a result of improper packaging of returned items
 - Lack of maintenance
 - Operational and normal wear and tear, which includes, for example, the following components:
 - Filters
 - Batteries
 - Articles for one-time usage, etc.
 - Failure to use original spare parts
- Weinmann is not liable for consequential harm caused by a defect if it is not based on intention or gross negligence. Weinmann is also not liable for minor physical injury to life or limb resulting from negligence.
- Weinmann reserves the right to decide whether to eliminate defects, to deliver a defect-free item or to reduce the purchase price by a reasonable amount.
- If Weinmann rejects a claim against the warranty, it assumes no expense for transport between customer and manufacturer.
- Implied warranty claims remain unaffected by these changes.

EEC Letter of Conformity on Medical Devices

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declare in sole
responsibility the

Product Name: **Home Mechanical Ventilator**

Type/Model: **VENTImotion**

being in conformity
with the respective
regulations of the
following guideline:

Medical Device Directive 93/42/EEC

Classification:

II a

Standards used

in particular:
EN 60601-1: 1990 + A1: 1993 + A2: 1995
EN 60601-1-2: 2001
DIN EN ISO 10651-6: 2004

Hamburg, 14.04.2005

Dr. K.-A. Feldhahn
Executive Vice President
Chief Operating Officer

01.2004 Alb

13. Index

Accessories	43
Alarms	35
Auto switch-on	9, 19, 20, 34
Carrying bag	23
Cleaning	25
Disinfection	25, 28
Dry hose system	26
Exhalation system	16, 19, 23, 25, 39, 44
Filter	34, 39, 45, 47
Function check	30
Inspiration time	9
Malfunctions	34
Mask	16, 19, 23, 25
Respiratory frequency	9, 47
Safety instructions	10, 12
Servicing	11
Side effects	13
Soft start	4, 19, 21, 22
Spare parts	45
Technical specifications	46
Therapy mode	9
Travelling	23
VENTI <i>click</i>	17, 22, 28, 46
VENTI-O ₂	17
VENTI <i>power</i>	18
Warranty	50

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Geräte für Medizin GmbH+Co.KG
Geräte für Medizin GmbH+Co. KG