ba12209de

pH 330/SET
pH 340/SET

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SenTix 41

pH-Einstabmeßkette
pH Combined Electrode

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Please read these safety guidelines attentively before the instrument is put into operation!

This instrument has been built and tested acc. to the German safety regulations for electronic test instruments IEC 1010 and has left our works in a condition complying with all the requirements of technical safety.

The perfect functioning and operational safety of the instrument is only ensured if the user observes the normal safety precautions as well as the specific safety guidelines stated in the present instruction manual.

- *Before connecting the plug-in power supply unit to the electricity supply network, ensure that the operating voltage stated on the plug-in power supply unit corresponds to the mains voltage (statement of the supply voltage range).*
- *The perfect functioning and operational safety of the instrument can only be maintained under the climatic conditions specified in the "Technical data" section of this instruction manual.*
- *When the instrument is moved from cold to warm surroundings, condensate may occur and interfere with the functioning of the instrument. In such a case, the user should wait until the temperature of the instrument has adapted to the ambient temperature before using the instrument again.*
- *Opening of the instrument, adjustment, maintenance or repair must only be carried out by a suitably qualified specialist authorized by us.*
- *If there is reason to assume that the instrument can no longer be employed without a risk, it must be set aside and appropriately marked to prevent further use.*
- *The safety of the user may be endangered, e. g. if the instrument*
 - *shows visible damage,*
 - *no longer operates as specified,*
 - *has been stored over a longer period under unsuitable conditions,*
 - *has been subjected to difficult conditions during transport.*
- *If in doubt, the instrument should as a rule be sent back to the manufacturer - "Wissenschaftlich-Technische-Werkstätten GmbH" for repair or maintenance.*

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SET equipment pH 330/340	10
Instruction manual pH 330/340	11
Control elements	12
Display	12
Keypad	14
Rear panel	14
Putting into operation	16
Battery operation (non-rechargeable battery versions only)	16
Exchange of batteries	16
Mains operation and operation with rechargeable batteries (rechargeable battery version only)	17
Display indicator of battery condition	18
Switching on / Display test	18
Calibration	20
Automatic calibration	20
Conventional two point calibration "ConCal"	22
Conventional single point calibration "ConCal"	25
Measurement	26
Redox voltage	27
Redox voltage	27
Check of the redox electrode ("Calibration")	27
Measurement	27
Temperature measurement	28
pH electrodes without temperature probe	28
Manual input of temperature value	28
Probe supervision	29
Probe symbol	29
Calibration interval	30
Drift control "AR"	31
Criteria	31
Calibration with drift control	31
Measurements with drift control	31
Data logging	33
Memory of measured values	33
Manual storing	33
Automatic storing	34
Memory full	36
Read-out of memory	37
Output of display	37
Erasure of memory	38


Analog output (Interface versions only)	39
Data transmission (Interface versions only)	39
Serial interface	40
Remote control	41
Calibration protocol	41
Automatic transmission of calibration protocols and measuring records	42
Timed Transmission of measuring records	42
Transmission of measuring records from the memory	43
Setup	45
RESET function	46
Troubleshooting	50
pH system messages	50
Other possible errors- pH measurement	51
General system messages	52
Accessories	55
General accessories	55
Accessories for pH measurement	56
Accessories for Redox measurement	56
Technical data	62
Test certificate EMC	65
Manufacturer's Certificate of Compliance: "Use in hazardous areas - zone 2"	66
Application reports	67
List of recommended application reports	67
 Instruction manual SenTix 41	
Operation / Maintenance	70
Putting into operation	70
Recommended areas of application	70
Cleaning	71
Storage	71
Replacement of electrolytic solution	71
Aging	72
Disposal	72
Trouble Shooting	73
Technical data	74

Instruction manual SenTix 21

Operation / Maintenance 76
 Putting into operation 76
 Recommended areas of application 76
 Cleaning 77
 Storage 77
 Replacement of electrolytic solution 77
 Aging 78
 Disposal 78
Trouble Shooting 79
Technical data 80

Instruction manual TFK 325/HC

Exchange of electrode 82
Recommended application range 82
Technical data 83

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- *KCl Storing solution for electrodes*
- *pH buffer solution STP 4, 50 ml*
- *pH buffer solution STP 7, 50 ml*
- *pH buffer solution STP 4, 50 ml*
- *Measuring beaker 50 ml*

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*pH combined electrode
SenTix 21 or SenTix 41*

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
*Line adaptor (rechargeable battery
version only)*

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pH pH value

User guidance and measured values:

pH, voltage, Slope of the pH electrode,
Asymmetry

Dimensions:

mV Voltage/Asymmetry
mV/pH Slope of the electrode

Calibration data:

Probe evaluation

Status:

TP Temperature measurement
active
°C Temperature


Measured values and set parameters:

Temperature, Time, Date, Numerator,
Number to identify measured value, Baud
rate

7	Status: RCL Funktion "Speicher Auslesen" aktiv	Status: RCL Funktion "Read-out of memory" active
8	Status: AR Driftkontrolle aktiv AR statisch Anzeige des stabilen Meßwertes AR blinkend Bestimmen des stabilen Meßwertes	Status: AR Drift control active AR static Stable measured value is displayed AR flashing Stable value is being determined
9	Kalibrierverfahren: AutoCal TEC AutoCal DIN ConCal	Calibration procedures: AutoCal TEC AutoCal DIN ConCal
10	Status: Store Funktion "Manuelles Speichern" aktiv Auto Store Funktion "Automatisches Speichern" aktiv	Status: Store Funktion "Manual storing" active Auto Store Funktion "Automatic storing" active
11	Status: LoBat Akkus/Batterien entladen	Status: LoBat Batteries discharged
12	Status: Time Uhrzeit Day, Month Tag und Monat Year Jahreszahl Baud Übertragungsgeschwindigkeit No. Nummer des Meßwert- speichers Ident Nummer zur Meßwert- identifikation	Status: Time Time Day, Month Day and month Year Year Baud Transmission speed No. Number of storage location Ident Number to identify measured value

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
Battery operation (non-rechargeable battery versions only)

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

Switch off the instrument before opening!

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Interior view:
Swing open lower shell

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Change all four batteries, pay regard to polarity!
(4 MIGNON, alkali-manganese batteries)

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*Use leakage-free alkali-manganese batteries only.
Take care that the gasket is not damaged when the instrument is closed.*

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Mains operation and operation with rechargeable batteries (rechargeable battery version only)

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*Use only rechargeable batteries
with the following specifications:
NiCd, AA size, 1.2 V/max.
750 mAh*

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*Operation time with one charging:
1500 hours.*

Charging time: 16 hours.

*The batteries are also charged when the
meter is switched off.*

*Normal measuring operation is also
possible during the charging.*

*Connect meter (socket 2) to mains via
line adaptor.*


**Line adaptor: Use only Friemann and
Wolf battery charger, model FW1199.**

	Input	Output
Friwo FW1199/11.7864, (Friwo Part. No. 17622613)	230 VAC (+10%/ -15%) 50 Hz / 5,6 VA	12 VDC / 130 mA / 1,56 VA
Friwo FW1199/11.7880, (Friwo Part. No. 1794043)	120 VAC (+10%/ -15%) 60 Hz / 6 VA	12 VDC / 150 mA

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
*Before the instrument is plugged in
and switched on, ensure that the
working voltage stated on the line
adaptor and the mains voltage
correspond.*

Batteriezustandsanzeige

Display indicator of battery condition

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When "LoBat" is displayed the batteries are nearly discharged. The meter is operable for no longer than another 10 hours.

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
The **economy circuit** switches off the meter approx. 1 hour after the last key operation.

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Exceptions:
Interface cable is connected or the function "Automatic storage" is switched on.

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

Stored value for the asymmetry e.g. 12 mV.

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After a total discharge set date and time if necessary.

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Admissible standard solutions:

WTW Technical Buffer Solutions
"Ed..."
Buffer solutions according to
DIN 19266 "Ed..." - "AutoCal":
A, C, D and F or
WTW models PL2, PL4, PL7 and
PL9.

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
Note:


Highest accuracy when the temperature
of standard solutions and sample solution
is the same.

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
Start measurement.
AR flashes.
Wait until the display shows:

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

– *Rinse pH electrode and immerse into the second Technical Buffer Solution.*

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
*Start measurement.
AR flashes.
Wait until AR extinguishes.*

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
*Offset voltage (Asymmetry).
Admissible range: ± 30 mV.*

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**Conventional two point
calibration "ConCal"**

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
Admissible standard solutions:

*Buffer solution pH 7.0 ± 0.5 and
any buffer solutions.*

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

*Immerse pH electrode into a
neutral buffer solution.
pH value of the solution:
7.0 ± 0.5*

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

Pay attention to the temperature of the buffer solution!

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

*Indication: Asymmetry.
Admissible asymmetry range:
 ± 30 mV*

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
Pay attention to the temperature of the buffer solution.

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
Indication: Slope of the electrode.
Admissible slope range:
-50 mV/pH ... - 62.0 mV/pH.

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
***Conventional single point
calibration "ConCal"***

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Admissible standard solution:
Standard solution suitable for the sample.

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

– Immerse pH electrode into the standard solution.

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– Set pH value of the standard solution, e.g. 1.49.

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
Pay attention to the temperature of the standard solution!

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
Indication: Asymmetry.
Admissible asymmetry range:
 ± 30 mV

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
Check of the redox electrode ("Calibration")

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
*Compare voltage and temperature value
to the data on the bottle of the redox buffer
solution.
If the difference is too high clean or
change electrode.*

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pH electrodes without temperature probe

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
For pH measurement it is also possible to use other electrodes than the SenTix 41. As most of these other electrodes have no or no suitable temperature probe the pH 330/340 offers the possibility to input the temperature manually for these applications:

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Set temperature value of sample solution.

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After a calibration the meter evaluates the condition of the calibrated pH electrode see also application report no. 497298: "Evaluation of electrochemical probes".

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Note:

S and U_{ASY} are evaluated separately. Then the worse result is used for the quality indication.

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
Action with evaluation "poor" or E3 see chapter "Troubleshooting"

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*When **probe symbol is flashing**:
Calibration interval has expired.
Calibrate pH electrode as soon as
possible.*

*The probe evaluation is made after the
calibration. The reliability of the meter,
therefore, depends on the calibration
interval considerably.
How to change the calibration interval is
described in the chapter "Setup".*

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*Calibration interval has expired.
Measurement still possible.
Calibrate as soon as possible.*

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Besides other factors, the stability of the measured value has a considerable influence on the reproducibility of the measured value.

The drift control checks the stability of the measuring signal and causes the stable measuring result to be displayed.

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Criteria

with unchanged measuring conditions

pH-Wert:	besser 0,02	pH value:	<i>better 0.02</i>
Einstellzeit:	> 30 Sekunden	Response time:	<i>> 30 seconds</i>

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Calibration with drift control

During calibration the meter automatically switches the drift control on or off respectively.


By pressing RUN/ENTER it is possible to cancel the drift control and to accept the value manually.

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

*Switch on drift control.
Measured value frozen.*

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
*AR display indicator flashes.
Wait until the measuring signal is
stable and AR stops flashing.*

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
*It is possible to cancel the drift
control and to accept the value
manually by pressing
RUN/ENTER.*

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Switch off drift control.

Meßwertspeicher

Memory of measured values

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Memory size: Maximum 200 measuring records

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
Storing measured data is only possible during measurement (measuring function switched on). The pH 330/340 always stores a complete measuring record consisting of:

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
Example: Storage of the 38th measuring record.

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*Last selected time interval.
Example 30 seconds.
With \blacktriangle or \blacktriangledown : Change the interval.*

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*Last selected identity number.
Input of a number to designate the
measured value (e. g. for the
sample location).*

*With \wedge or \vee : change the identity
number.*

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*Indication: Number of storage
locations that are still
free.*


*Automatic storage is started.
All functions are blocked.*

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
*After the automatic storing has
occupied all the 200 storage
locations the meter switches off
automatically.*

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
*After occupation of the 200th
storage location the meter displays
Sto Full.*

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
Keep corresponding key depressed:
Scrolling.

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One seconds no key pressed:
*Stored temperature is indicated instead of
the search criterion. To go on press
corresponding key.*

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
*The erasure function erases **all the stored measuring records**. Calibration protocols (Interface versions only) remain stored up to the next calibration.*

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*Erasure function switched on.
(Only appears when the memory contains data.)*

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*Confirm erasure process.
Any other key cancels the erasure process.*

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
Analog output
(Interface versions only)

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Connect analog output to recorder with cable AK 323. The output is automatically switched over to recorder output.

Connections:

- 1 free
 - 2 plug codification
 - 3 ground
 - 4 analog output
- Connection of 2 and 3 activates analog output.*

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The connected recorder cable AK 323 blocks the automatic switch-off.

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Data transmission
(Interface versions only)

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The output of calibration protocols and measuring records to the printer or computer is made via the serial interface.

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
Connect serial output of the meter to serial connection of the PC/printer.

PC: Cable AK 340/B,
Order no. 902 841

Printer: Cable AK 325/S,
Order no. 902 837

The input/output is automatically switched over to RS 232.

Connections:

- 1 CTS
- 2 Plug codification + RxD
- 3 Mass
- 4 TxD

Connecting 2 and 3 and connection of a PC/printer activates the serial output.

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Set transmission data at the PC/printer:

Baud rate: 1200, 2400,
4800, 9600 selectable
Parität (PC only): none
Handshake: RTS/CTS + Xon/Xoff
Data bits (PC only): 8
Stop bit (PC only): 1

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Economy circuit:

The meter automatically switches itself off approx. 1 h after the last key operation.

Exception:

Continuous operation on PC or printer with connected interface cable and activated timed output.

For setup of the serial interface (Baud rate, timed transmission, transmission on keypressing) please refer to chapter "Setup".

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Remote control

The meter can be operated by remote control from a PC via the AK 340/B cable.

The meter is operated with commands simulating key operations or requesting the current display contents. A more detailed description of the individual remote control commands and example programs can be found in the communication set KOM 340 (order no. 902 851).

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Calibration protocol

- *Printing date Printing time*
- *Series number of the printed board (to identify the used measuring electronics)*
- *Heading "CALIBRATION pH"*
- *Date / time of last calibration*
- *Calibration interval*
- *Used calibration procedure/ temperature auto/man.*
- *Used standard solutions (marked with *)*
- *Voltage 1, Temperature 1*
- *Voltage 2, Temperature 2*
- *Slope value S*
- *Offset voltage ("Asymmetry": ASY)*
- *Probe evaluation*

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```
CALIBRATION PROTOCOL
21.03.97            11:25
Device No.: 02307424
CALIBRATION pH
Cal Time: 21.03.97 / 11:25
Cal Interval: 10 d
AutoCal TEC      Tauto AR
Buffer 1           2.00 *
Buffer 2           4.01
Buffer 3           7.00 *
Buffer 4           10.00
C1                13 mV    22.3 °C
C2                298 mV   21.9 °C
S                 -58.7 mV/pH
ASY               3 mV
Sensor            +++
```

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**Automatic transmission of
calibration protocols and
measuring records**

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After calibration the meter automatically transmits the calibration protocol or, after measurement with drift control (AR), the corresponding measuring record to a connected printer or PC via the serial interface (see chapter "Calibration protocol").

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
**Timed Transmission of measuring
records**

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
*When **Int.2 = OFF** is set a **single** measuring record is transmitted.*

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*Last selected time interval.
Example 30 seconds.
With \blacktriangle or \blacktriangledown change the interval.*

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Start of the transmission with the selected time interval.

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Example for timed output during measuring operation:

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Transmission of measuring records from the memory

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Each transmission starts with the output of the calibration protocol.

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

Start of the transmission of the calibration protocol.
Termination with RUN/ENTER or pH/mV.

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The transmission of the measuring records starts with the oldest value.

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*The output of the stored data is also possible if **hlt 2 OFF** is set (but the cable must be connected).*

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The output can be terminated any time (e.g. after printing the calibration protocol)

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No. 3:
12.5.97 14:24
+ 1210 mV 38.9 °C
Tauto
Ident: 87

No. 4:
13.5.97 8:17
pH 7.03 18 °C
Tman AR
Ident: 55

*Record number 1
Date/time of storage
pH value with temperature
Autom. temperature
Number of sample location*

*Record number 2
Date/time of storage
pH value with temperature
Autom. temp., AutoRead
Number of sample location*

*Record number 3
Date/time of storage
mV value with temperature
Autom. temperature
Number of sample location*

*Record number 4
Date/time of storage
pH value/temperature
Manual temp., AutoRead
Number of sample location*

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
Changing the following parameters is possible in the sequence listed below (Default setting = printed in bold):

Anzeige / Display	Parameter / Parameters	Einstellmöglichkeiten / Setting possibilities
Baud	Baudrate (nur Schnittstellenversionen)	Baud rate (Interface versions only)
Time Int. 1	Intervallzeit (Speichern)	Interval times (storage)
Time Int. 2	Intervallzeit ser. Übertragung (nur Schnittstellenversionen)	Interval times ser. transmission (Interface versions only)
pH Int. 3	Kal.-Intervall pH	Cal. interval pH
Day/Month	Datum (Tag)	Date (day)
Day/Month	Datum (Monat)	Date (month)
Year	Datum (Jahr)	Date (year)
Time	Uhrzeit (Stunde)	Time (hour)
Time	Uhrzeit (Minute)	Time (minute)

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

*Resetting of single measuring functions to default condition.
Even after a RESET the following data remain stored:*

- *Measured values stored in memory*
- *Date/Time*
- *Calibration interval Int.3*

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For each measuring function select with or :

- Yes: Perform RESET
No: Do not perform RESET*

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
pH system messages

System message	Cause	Action
<p>QFL Indication range exceeded</p>	<p>pH electrode:</p> <ul style="list-style-type: none"> • Not immersed • Cable broken • Electrolytic gel dried • Air bubble in front of diaphragm • Air in diaphragm 	<p>pH electrode:</p> <p>Immerse into sample solution</p> <p>Exchange electrode</p> <p>Exchange electrode</p> <p>Remove air bubble</p> <p>Suck off air or moisten diaphragm</p>
<p>E3 Slope or asymmetry outside of admissible ranges (Slope: -50 ... -62 mV/pH; Asymmetry: -30 ... +30 mV)</p>	<p>pH electrode:</p> <ul style="list-style-type: none"> • pH electrode too old • Diaphragm contaminated • Membrane contaminated • Moisture in the plug • Electrolytic gel too old <p>Meter:</p> <ul style="list-style-type: none"> • Wrong calibration procedure • Wrong temperature of solution (Without temperature probe only) • Moisture in the socket <p>Buffer solution:</p> <ul style="list-style-type: none"> • Unsuitable buffer solutions • Buffer solutions too old • Used up buffer solutions 	<p>pH electrode:</p> <p>Exchange pH electrode</p> <p>Clean diaphragm</p> <p>Clean membrane</p> <p>Dry plug</p> <p>Exchange electrode</p> <p>Meter:</p> <p>Select correct procedure</p> <p>Set correct temperature</p> <p>Dry socket</p> <p>Buffer solution:</p> <p>Use WTW Technical buffer solutions</p> <p>Pay attention to perishability of solutions.</p> <p>Use solutions only once.</p>

Other possible errors - pH measurement

System message	Cause	Action
<i>Voltage of pH electrode not stable</i>	<p>pH electrode:</p> <ul style="list-style-type: none"> • Diaphragm contaminated • Membrane contaminated <p>Sample solution:</p> <ul style="list-style-type: none"> • pH value not stable • Temperature not stable <p>Electrode + sample solution:</p> <ul style="list-style-type: none"> • Conductivity too low • High temperature • Organic fluids 	<p>pH electrode:</p> <p>Clean diaphragm Clean membrane</p> <p>Sample solution:</p> <p>Measure in an airproof vessel if necessary Temper if necessary</p> <p>Electrode + sample solution:</p> <p>Use suitable pH electrode Use suitable pH electrode Use suitable pH electrode</p>
<i>Measured values obviously wrong</i>	<p>pH electrode:</p> <ul style="list-style-type: none"> • Unsuitable pH electrode • Temperature difference between calibration and sample solution too high • Interference voltages <p>• Unsuitable measuring procedure</p>	<p>pH electrode:</p> <p>Select suitable pH electrode Temper buffer or sample solution</p> <p>With grounded sample solution do not use a grounded PC or printer</p> <p>Pay attention to special procedures</p>

General system messages

System message	Cause	Action
<i>Meter does not react on keypressing</i>	<i>Undefined operation state of the processor, e. g. after inadmissible EMC admission (outside the specified ranges)</i>	<i>Reset of processor: Press RCL key and switch meter on.</i>
LoBat	<i>Batteries nearly discharged.</i>	<i>Charge rechargeable batteries/ Change non-rechargeable batteries.</i>
	<i>Timeout of serial interface.</i>	<i>Check connected instrument.</i>
 <i>is flashing</i>	<i>Set calibration interval has run out.</i>	<i>Recalibrate probe.</i>

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General accessories

	Model	Order no.
Adaptor cable <i>for connection to a recorder via analog output</i>	AK 323	902 835
Interface cable <i>for connection to a printer via RS 232 interface</i>	AK 325/S	902 837
Interface cable <i>for connection to a PC via RS 232 interface, incl. adaptor 9/25 pin and operation manual for current terminal programs</i>	AK 340/B	902 841
Communication set <i>consisting of Interface cable AK 340/B, with adaptor 9/25 pin and transmission software Multi/ACHAT II</i>	KOM 340	902 851
Telescop <i>1.5 - 3 m with sensor holder for pH combined electrode SenTix...</i>	TH 340	402 802
Protection armoring <i>Shock absorbent, rubber elastic, with stand and clamping holder for sensor cable</i>	SM 325	902 818
Beaker set <i>for protection armoring SM 325, consisting of one electrode beaker pH, holding device and carrying strap for field application</i>	TG/pH	902 830
Field armoring <i>Rubber elastic protection armoring with integrated carrying handle, electrode beaker pH, carrying strap, including holding device and stand for measurements in the laboratory</i>	FM 325/pH	902 832

Accessories for pH measurement

	Model	Order no.
pH combined electrode with integrated temperature probe and watertight plug (IP 67)		
with cable lengths 1 m	SenTix 41	103 635
3 m	SenTix 41-3	103 636
pH combined electrode with watertight plug (IP 67)		
with cable length 1 m	SenTix 21	103 631
3 m	SenTix 21-3	103 632
Holding clip with integrated temperature probe		
for SenTix 21	TFK 325/HC	109 818
for SenTix 21-3	TFK 325/HC-4	109 819
Technical buffer solutions, 250 ml		
Buffer pH 4.01	TPL 4	108 800
Buffer pH 7.00	TPL 7	108 802
Buffer pH 10.00	TPL 10	108 804
DIN/NBS buffer solutions 250 ml		
buffer solution pH 1.679	PL 2	109 000
buffer solution pH 4.006	PL 4	109 110
buffer solution pH 6.865	PL 7	109 120
buffer solution pH 9.180	PL 9	109 130
Storing solution		
250 ml KCl 3 mol/l	KCl-250	109 705

Accessories for Redox measurement

	Model	Order no.
Redox combined electrode (platinum - silver/silver chloride), immersion length 120 mm, with plug head	SenTix ORP	103 648
Connection cable for electrodes with DIN plug, 1 m cable	AS/DIN	108 110
Redox buffer solution for control of Redox electrodes $U_{H^+} = 427$ mV, 1 bottle with 250 ml	RH 28	109 740
Reagents for regeneration of Redox electrodes consisting of 10 g activation powder and 30 g chlorina powder	SORT / RH	109 730

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Housing

Length	172 mm
Width	80 mm
Height	37 mm
Weight	approx. 0.3 kg (total weight of meter)
Color	black, RAL 9005

EMC and VDE Standards

Protection degree	IP 66	IEC 529 (without line adaptor)
Protection class	3	IEC 1010
Climate class	JVD	DIN 40040
EMC Emission	Generic standard EN 50081-1 standard part FCC class A	
EMC Immunity	Generic standard EN 50082-1 standard part	

Operation conditions

Operating temperature	-10 ... 55 °C
Storing temperature	-25 ... 65 °C
	Operation in the open, dew possible
Relative humidity	Yearly mean: < 90 % 30 days/year: <100 % Other days: < 80 %

Power supply

Non-rechargeable batteries (non-rechargeable battery versions only)	4 x 1.5 V batteries typ AA or
Rechargeable batteries (rechargeable battery version only)	Rechargeable batteries: 4 x NiCad Typ AA, 1.2 V/max. 750 mAh
Working time	1500 operation hours with NiCad rechargeable batteries 5000 operation hours when using alkali/ manganese batteries
Mains (rechargeable battery version only)	Line adaptor FRIWO FW 1199/11.7864, Friwo Part No. 1762613 Input: 230 VAC/(+10%/-15%)/50Hz/5.6 VA Output: 12 VDC/130 mA/1.56 VA Line adaptor FRIWO FW 1199/11.7880, Friwo Part No. 1794043 Input: 120 VAC/(+10%/-15%)/60Hz/6 VA Output: 12 VDC/150 mA Line adaptor with UK plug

Data security

unlimited, not depending on charging condition of the batteries

Time

Clock	Real time clock
Accuracy	1 min/month max. deviation

Memory of measured values	
Type	Ring store
Quantity	200 measuring records
Timer for timed storing	5 sec, 30 sec, 1 min, 5 min, 10 min, 15 min, 30 min, 60 min
Outputs (Interface versions only)	
The interface (RS 232 or recorder) is automatically recognized and switched over accordingly when the meter is connected.	
Analog output	Automatic switch-over when the recorder cable AK 323 is connected. Output signal corresponds to electrode voltage, i.e. 50 - 62 mV/pH unit Accuracy: $\pm 0.1\%$ of meas. value ± 4 mV at ambient temperature $-10^{\circ}\text{C} \dots +55^{\circ}\text{C}$ Internal resistance: < 5 Ohm (current limit max. 0.2 mA output current)
Serial interface (Interface versions only)	Automatic switch-over when the cable AK 340/B or AK 325/S is connected.
Type	RS 232, data output
Baud rate	adjustable 1200, 2400, 4800, 9600 Baud
Data bits	8
Stop bit	1
Parity	None
Handshake	RTS/CTS + Xon/Xoff
Cable length	max. 15 m
Output format	ASCII strings (IBM font) Field separation by CR, LF
Timer for timed output	5 sec, 30 sec, 1 min, 10 min, 15 min, 30 min, 60 min, OFF
Connections (electromechanical construction)	
pH/voltage	Socket according to DIN 19262 (Protection degree IP 67) in connection with SenTix electrode
Temperature	1 pin WTW jack (Protection degree IP 67)
Line adaptor (rechargeable battery version only)	2 pin socket with pin contacts for FRIWO plug
Serial interface, analog output (Interface versions only)	4 pin socket with pin contacts

pH measurement	
<i>(Pay attention to the technical data of the electrode additionally)</i>	
Measuring range	-2.00 ... 16.00
Resolution	0.01
Accuracy	0.01 ± 1 digit
General points	Input resistance: > 5 · 10 ¹² Ohm Input currence: 5 · 10 ⁻¹³ A typically < 2 · 10 ⁻¹² A at 25 °C
Usable electrodes	pH electrodes with plug according to DIN 19262 and electrode zero point pH = 7.0 ± 0.5
Calibration procedure	<u>two-point-calibration:</u> AutoCal TEC WTW Technical Buffer Solutions AutoCal Buffer Solutions according to DIN 19 266 <u>one-point-calibration:</u> ConCal 1 or 2 point calibration with suitable buffer solutions
Calibration interval	1 ... 999 days (preset 7 days)
Slope range	-62 ... -50 mV/pH
Asymmetry	-30 ... +30 mV
Redox measurement	
<i>(Pay attention to the technical data of the electrode additionally)</i>	
Range	-1250 mV ... + 1250 mV
Resolution	1 mV
Accuracy	1 mV ± 1 digit
Usable electrodes	Pt-Electrodes with plug according to DIN 19262
Temperature	
<i>(Pay attention to the technical data of the probes additionally)</i>	
Measuring range	-5 ... 99.9 °C
Resolution	0.1 K
Accuracy (meter)	with NTC 30: ± 0.1 K ± 1 digit at operating temperature -10 °C ... + 55 °C with Pt 1000: <u>accuracy at operating temperature</u> -10 °C ... +15 °C ± 0.5 K ± 1 digit 15 °C ... 35 °C ± 0.1 K ± 1 digit 35 °C ... 55 °C ± 1 K ± 1 digit
Usable probes	30 KOhm NTC integrated in the electrodes; TFK 150 or TFK 325 as separate probe; Pt 1000 as separate probe or integrated in electrodes.
man. Temperature input	-20 ... 130 °C when no temperature probe is connected.

Test certificate EMC

SENTON

EMV-Prüfzentrum EMI/EMC-Testcenter

CERTIFICATE
to the Electromagnetic Compatibility

to Test Report No. 52501-70351-1/-2

Sample:	pH 3XX pH meter
Device:	pH 340-A
Uniform EMC design:	pH 330, pH 340, pH 340/ION
Applicant:	WTW GmbH
Regulations:	EN 50081-1:1992 EN 50082-1:1992 FCC part 15 subpart B

Testresult:

The samples are in compliance with the RFI requirements and the immunity requirements according to above referenced regulations. The following severity levels have been achieved:

RFI Emissions

Requirements according to EN 50081-1:1992
Requirements according to FCC part 15 subpart B limit class A

Immunity

Electrostatic Discharge IEC 801-2:1984	Air discharge	8 kV
Electromagnetic Fields IEC 801-3:1984	27 MHz - 500 MHz	3 V/m
Electrical Fast Transients IEC 801-4:1988	AC power input port DC and data/signal lines	1 kV 0.5 kV

Straubing, May 30, 1997

Senton GmbH
Johann Roidt

**Manufacturer's Certificate of Compliance:
"Use in hazardous areas - zone 2"**

Manufacturer's Certificate

We hereby confirm that the

**pH/mV Meter
pH 330, pH 340**
(without power-supply unit)

*can be utilized in accordance with
the regulations DIN/VDE 0165
in explosion endangered areas of zone 2.*

*This concerns the equipment where
no sparks, electric arcs or undue temperatures
occur at operative use.*

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List of recommended application reports


	Appl. report no.
pH value	
<i>pH measurement in wastewater</i>	1296295
<i>Influence of temperature on pH measurement</i>	1095244
<i>Stability of pH buffer solutions</i>	196065
Redox voltage	
<i>Redox measurement in wastewater</i>	1296297

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pH combined electrode with integrated temperature probe

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Some of the potassium chloride solution may leak out of the wetting cap during transport or storage, leaving a crust of potassium chloride solution after drying. This layer of salt is harmless and can be rinsed off with water.

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In-situ measurement and simple laboratory measurement in aqueous media. The durability can be shortened if using the combined pH electrode near the limits or outside the specified areas of application (see technical data).

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<i>Impurities</i>	<i>Cleaning agents</i>	<i>Contact time</i>
<i>Water-soluble impurities</i>	<i>Distilled water</i>	<i>10 minutes</i>
<i>Other impurities</i>	<i>Warm water with detergent</i>	<i>10 minutes</i>

Subsequently, thoroughly rinse with deionized water.

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Aging

Please note that each pH electrode undergoes a natural aging process, as a result of which the reading will become increasingly slow while the slope of the pH electrode gradually decreases.

Incorrect storage, certain measurement media (e.g., sulphate-containing, highly alkaline, protein-containing media) or high temperatures as well as large pH or temperature changes will shorten the service life considerably.

The guarantee does not apply in the case of failure caused by the measurement media or mechanical damage.

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Disposal

Your old pH electrodes can be returned to WTW who will ensure correct disposal in compliance with the legal regulations.

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
<i>Process</i>	<i>Result</i>	<i>Action</i>
<i>Calibrating</i>	<i>Error message (calibration not possible)</i>	<ul style="list-style-type: none"> ● <i>Remove wetting cap</i> ● <i>Clean</i> ● <i>Check plug-in contact of DIN plug</i> ● <i>Check buffer solutions</i>
<i>pH measurement in a solution</i>	<i>Error message (overrange) Drifting or unstable values</i>	<ul style="list-style-type: none"> ● <i>Check plug-in contact of DIN plug</i> ● <i>Clean</i> ● <i>Dry plug-in contact of DIN plug</i> ● <i>Only store pH electrode in the specified solution</i> ● <i>Conductivity of measurement solution too low or changing pH value</i>
<i>Temperature measurement</i>	<i>Error message (Exceeding of meas. range) Drifting values</i>	<ul style="list-style-type: none"> ● <i>Control plug contacts</i> ● <i>Dry plug contacts</i> ● <i>Wait for temperature conformity</i>

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pH-Einstabmeßkette
pH-combined electrode

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Putting into operation

- *Remove wetting cap.*
- *Soak in neutral buffer solution for approx. 10 minutes.*
- *Calibrate according to instrument instructions.*

If the membrane has dried out, soak the pH electrode in a neutral buffer solution or, preferably, potassium chloride solution ($c = 3 \text{ mol/l}$) for 24 hours.

Some of the potassium chloride solution may leak out of the wetting cap during transport, leaving a crust of potassium chloride solution after drying. This layer of salt is harmless and can be rinsed off with water.

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Recommended areas of application

In-situ measurement and simple laboratory measurement in aqueous media. The durability can be shortend if using the combined pH electrode near the limits or outside the specified areas of application (see technical data).

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Aging

Please note that each pH electrode undergoes a natural aging process, as a result of which the reading will become increasingly slow while the slope of the pH electrode gradually decreases. Incorrect storage, certain measurement media (e.g., sulphate-containing, highly alkaline, protein-containing media) or high temperatures as well as large pH or temperature changes will shorten the service life considerably. The guarantee does not apply in the case of failure caused by the measurement media or mechanical damage.

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
Process	Result	Action
<i>Calibrating with buffer solutions pH = 7 and pH = 4</i>	<i>Error message (calibration not possible)</i>	<ul style="list-style-type: none"> ● <i>Clean</i> ● <i>Remove wetting cap</i> ● <i>Check plug-in contact of DIN plug</i> ● <i>Check buffer solutions</i>
<i>pH measurement in a solution</i>	<i>Error message (overrange) Drifting values</i>	<ul style="list-style-type: none"> ● <i>Check plug-in contact of DIN plug</i> ● <i>Clean</i> ● <i>Dry plug-in contact of DIN plug</i> ● <i>Only store pH electrode in the specified solution</i> ● <i>Conductivity of measurement solution too low or changing pH value</i>

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In-situ and laboratory measurements in drinking-, surface-, ground- and wastewater.

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