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MULTIMAT[®] C

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso

DENTSPLY

DENSPLY

MULTIMAT® C

* ↓ ↑ ↺ ↻

50 0500 max 0940 t0430
01 01 U 00.5 01.0 0

⊙ ☹ ☹ ☹ ☹ ☹ ☹ *

on/off	M	start/stop	7	8	9	R
▲			4	5	6	↑
◀		▶	1	2	3	↓
▼			C	0	E	S

MULTIMAT® C

MULTIMAT® C

D	Gebrauchsanweisung	7–23
GB	Instructions for use	25–41
F	Mode d'emploi	43–59
I	Istruzioni per l'uso	61–77
E	Instrucciones para el uso	79–95

Inhalt

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Contents

1 List of parts	27
1.1 List of parts	27
1.2 Technical data	28
2 Introduction and performance data	28
2.1 Introduction	28
2.2 Performance data	28
3 Safe use	28
3.1 Intended use	28
3.2 Hazards and safety notes	28
4 Setup and first operation	29
4.1 Unpacking	29
4.2 Checking accessories	29
4.3 Setup	29
4.4 First operation	29
4.5 Test program	29
5 Practical application: introduction	30
5.1 Limit values of the firing parameters	30
5.2 Display	30
5.3 Function	31
5.4 Keypad functions	31
5.5 Menues	32
6 Create individual programs	33
6.1 Load program	33
6.2 Enter preheating temperature	33
6.3 Enter heating rate	33
6.4 Enter firing temperature	33
6.5 Enter drying time	33
6.6 Enter preheating time	33
6.7 Enter vacuum level	33
6.8 Enter vacuum time	33
6.9 Enter firing time	34
6.10 Save program	34
6.11 Start program	34
6.12 Copy program	34
6.13 Special program	34

7 Specialised functions	34
7.1 Standby	34
7.2 Vacuum	34
7.3 Rapid cooling	34
7.4 View function	34
7.5 Power failure backup	35
7.6 Testing the firing chamber temperature with the silver test	35
7.7 Connecting a printer	35
8 Servicing and maintenance	36
8.1 Lift	36
8.2 Firing plynth with firing plynth holder	36
8.3 Vacuum pump	36
8.4 Replacing the firing muffle	36
8.5 Replacing the control element	37
9 Faults and remedies	37
10 Preset programs	39
11 Warranty	40
12 EU Declaration of Conformity	41

Information symbols

Warning!

This symbol indicates important notes and instructions; non-compliance could result in injury and accident hazard.


Caution!



This symbol indicates important information and notes that you should follow.

1 List of parts

1.1 List of parts; appliance view.

- 1. Cover
- 2. Housing
- 3. Firing plinth
- 4. Lift
- 5. Workpiece shelf
- 6. Control
- 7. Mains fuse
- 8. Connector for vacuum tube
- 9. Power connection
- 10. Electrical connection for vacuum pump
- 11. Transformer fuse
- 12. Firing muffle
- 13. Vacuum pump

1.2 Technical data

Diameter:	280 mm
Height:	590 mm
Voltage:	100, 110, 115, 125, 230 Volt
Frequency:	50/60 Hz
Performance:	1300 VA without pump
Weight:	17.5 kg

2 Introduction and performance data

2.1 Introduction

Dear customer,

Thank you for purchasing the Dentsply Multimat C.

This furnace represents an advanced product of the Multimat series.

Simple and direct operation allows for quick training times.

The Multimat C can store up to 99 customised programs.

7 factory-preset programs can be allocated to program numbers 92 – 98 if required.

All firing data is shown on an illuminated display.

All firing processes are controlled by microprocessors.

Due to its small base, the Multimat C can be set up even in limited space.

The instrument complies with all valid EU directives and VDE safety regulations.

Please read the operating manual !

2.2 Performance data

- Quartz muffle with freely radiating heat coil. Useable diameter 93 mm, height 67 mm.
- 2 x 24-digit LCD display capable of displaying alphanumeric characters
- Complete display of firing data
- Membrane keypad
- 99 fully customizable programs
- 7 factory-preset programs (Ceramo 3)
- Intervene to alter a running program
- Change a saved program to an individual special program
- Copy a saved program to another program number
- Firing temperature up to 1200 °C
- Vacuum time up to 99.5 minutes
- Display for theoretical and actual vacuum level
- Adjustable vacuum
- Firing time up to 99.9 minutes
- Drying and preheating time up to 25 minutes
- Temperature increase rate 1 – 120 °C/minute
- View function
- 3 cooling stages
- Rapid cooling programmable and manual
- Display of remaining time until end of firing
- Overheating protection

- Temperature display in Fahrenheit/Celsius
- Muffle monitoring
- High temperature accuracy
- Standby switch for preventing moisture accumulation in the firing chamber
- Display of hours of operation
- Display of vacuum pump running time
- Language selection (5)
- Error messages
- Automatic continuation of firing after brief power failure
- Unlimited data backup during power failure
- Automatic frequency change
- Acoustic signal, can be disabled
- Correction factor can be entered during the silver test
- Vacuum program to prevent moisture in the fibre insulation
- GS and CE marks

3 Safe use

3.1 Use within specifications

Caution !

The Multimat C was developed and is intended solely for the purpose of firing dental ceramics.

Dentsply is not liable for any damage resulting from non-specified use.

Intended use also includes all notes in this operating manual as well as all notes in the separate operating manual of the manufacturer of the vacuum pump.

3.2 Hazards and safety notes

Caution !

In order to ensure safe operation, please observe the following notes:

Do not set up furnace and vacuum pump next to heat sources.

The distance to the next wall should be at least 25 – 30 cm.

The surface should be non-flammable, and there should be no flammable objects nearby.

Place the vacuum pump in a well ventilated place. Vacuum pumps lubricated with oil should always be set up in a lower place than the furnace.

Do not touch parts that become hot, particularly the cover.

Prior to powering up, ensure that the mains power corresponds to the voltage specified on the instrument.

After any extended inactivity period, at high humidity or low temperatures, it may no longer be possible to generate a sufficient vacuum. In this case, preheat the furnace to approx. 1000 °C and afterwards allow to cool down to idle temperature with opened firing chamber using the rapid cooling function.

If the furnace is under vacuum in "Off" state for an extended period of time, the O-ring of the lift plate may adhere slightly.

When heating the firing muffle, there may be vibration noises from the heater coil.

Warning!

The instrument may only be connected to a **grounded socket** with a < 30 mA **leakage current** circuit breaker.

Always disconnect the instrument from main power when carrying out repairs and maintenance work.

Repairs on opened unit may only be carried out by **trained personnel**.

After repairs, a high voltage and protective ground wire test must be carried out.

In the case of faults or damages that prevent safe operation, the instrument must be secured against unintentional use.

4 Setup and first operation

4.1 Unpacking

- Check the package and read the unpacking instructions lying on top before proceeding.
- Unpack furnace and accessories according to the unpacking instructions.

4.2 Check accessories

- Check if the scope of delivery is complete.
- Each Multimat C is supplied with:
 - 1 workpiece carrier W incl. pins
 - 1 workpiece shelf
 - 1 firing plynth
- Inspect delivery for possible transport damage.
- If the delivery is incomplete or damaged during transport, please contact your supplier.

4.3 Setup

- Place the furnace on a suitable surface and allow for sufficient distance to walls (min. 25 cm).
- Connect the power plug of the Dentsply vacuum pump with the vacuum pump socket and insert the vacuum tube in the tube connector. The arrows on the filter must point into the direction of the vacuum pump.

4.4 First operation

Before completing the connection, check whether the mains voltage of your power supply corresponds to the mains voltage specified on the type label.

If the voltage does not correspond, the furnace must not be connected.

- Connect the power lead of the furnace with a mains socket protected by a 16 A fuse and a leakage current circuit breaker (< 30 mA).

- The firing chamber will now be ventilated and the language selection appears on the screen:

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- **Please choose your language within one minute by pressing the appropriate number.**
- The furnace now runs through a short initialization phase. After completion, it displays the version number of the software for 3 seconds. An automatic memory test follows and either program 99 or the last set program is loaded as test program.
- **Please press the button Lift ↓ immediately, the lift moves to the lower end position. Afterwards place firing plynth on its holder.**

Notes:

If no language is selected within one minute, the furnace automatically selects the language last used, generates a vacuum of approx. 800 hPa, if set in the menu, and switches to "Off".

Display (switched "OFF")

8	OFF	8
----------	------------	----------

(Theoretical value 800 hPa)

(Actual value approx. 800 hPa)

When in "Off" state, the furnace can be switched on with the button "on/off". The firing chamber will be ventilated and the program last used is displayed.

4.5 Test program

This test program will provide you with a first impression of the functions of the Multimat C. In factory condition, this program is saved under program number 99 and contains the following data:

Program number	99
Theoretical value of preheating temperature:	500 °C
Temperature increase rate	80 °C/minute
Theoretical value of firing temperature	940 °C
Temperature mode	c = Celsius
Actual value of firing temperature	Displays the actual temperature in the firing chamber
Drying	1 minute
Preheating	1 minute
Vacuum level	50 hPa
Vacuum time	0.5 minutes
Firing time	1 minute

The data of the test program can be changed.

Note: When heating the firing chamber from room temperature, there may be vibration noises from the heater coil for some seconds.

They are generated from the free movement of the heater coil in the quartz tube and are not significant.

5 Practical application: introduction

5.1 Limit values of the firing temperature



Parameter	Lower limit	Upper limit
Preheating temperature - theoretical	100 °C	1200 °C
Heating rate	1 °C/min.	120 °C/min. or max.
Firing temperature - theoretical	300 °C	1200 °C
Firing temperature - actual	0 °C	1250 °C
Drying time	0 min.	25 min.
Preheating time	0 min.	25 min.
Vacuum level	1 hPa	99 hPa
Vacuum time	0.00 min.	99.5 min.
Firing time	0.00 min.	99.5 min.
View	00.1 Min	3.0 min.
Standby temperature	120 °C fixed	
Standby vacuum	800 hPa fixed	
Cooling	3 stages	(1, 2 and 3)

Values outside these limits are not accepted when saving, starting and altering a program during firing. A signal will sound to inform the user.

5.2 Display

Illuminated display for simultaneous display of all firing parameters.

Key to the symbols and the appropriate firing temperature:

Start test program

Press button "start/stop".

The start control light illuminates and a blinking arrow, left of the drying time, indicates that the preheating time has not been reached.

After reaching the preheating time, the approximate total firing time/remaining firing time is displayed in the lower left corner and the program sequence starts with the drying phase.

A flashing of the data display of the program step indicates that this step is being processed and the appropriate time is counted down to zero.

At the end of every program step the relevant data display goes out. The remaining total firing time/remaining firing time is newly approximated after each program step has been completed.

Once the set vacuum level is reached, the temperature increase starts.

When the test program is finished, the lift moves down and a acoustic signal sounds.

5.3 Function

Program number

Shows the current program number. Up to 99 individual programs can be stored.

Rapid cooling (programmed)

Rapid cooling can also be saved in the menu together with the program. The symbol “❄” is then displayed under the symbol with the three ice crystals.

During rapid cooling, the vacuum pump switches on and aspirates air through the firing chamber until the starting temperature has been reached again.

Rapid cooling (manual)

Switch on by pressing the button Lift “↓”.

The function can only be switched on if:

- the program has ended
- the firing chamber is completely opened
- the firing chamber temperature is higher than the theoretical preheating temperature.

Preheating temperature - theoretical

100 – 1200 °C

Heating rate

1 – 120 °C/minute

Max. = full power (approx. 140 °C)

Firing temperature - theoretical

300 – 1200 °C

Actual temperature

0 – 1250 °C

Total firing time / error display

Shows the total firing time as approximate addition of all program steps. During the program sequence, the total firing time is updated after each completed program step. This also serves as the error display.

In this case, the error message overwrites the total firing time.

Drying time

0 – 25 minutes

Preheating time

0 – 25 minutes

Vacuum display

The vacuum display activates as soon as the cursor is placed on the display. Vacuum level 1 – 99 hPa.

09–01 flashing = vacuum is being generated 900 – 100 hPa.

< = 99 hPa flashing – absolute vacuum

Once the desired vacuum level is reached, the temperature increase to the set firing temperature starts.

Vacuum time

0 – 99.5 minutes

Firing time

0 – 99.9 minutes

Cooling stages

0 = Lift moves immediately to lower end position – no cooling

1 = Workpiece moves approx. 7 cm from the firing chamber

2 = Workpiece moves approx. 5 cm from the firing chamber


3 = Workpiece remains in the firing chamber

If programmed, the cooling starts after the end of the firing time.

It lasts until the preheating temperature (with firing muffle switched off) is reached again.

The 3 cooling stages of the Multimat C reduce tension within the ceramic.

5.4 Keypad functions


Button “on/off”

By pressing this button, the furnace control is switched on and off.

In “Off” state, the firing chamber is (will be) closed. Afterwards, if set in the menu “Configuration”, a vacuum of 800 hPa is generated or the standby temperature is activated (the pump must remain operational), the heating muffle switches off and the display shows the “Off” state, as illustrated below. This prevents moisture from entering the insulation.

Display:

Menu

Press this button to call up the menu selection.

Button start/stop

Press this button to start or abort a program.

Control light:

Illuminates while a program is running.

Caution !

The command "Start" is accepted in any lift position. If the steps "Drying" and "Preheating" are selected, the order to start is saved until the preheating temperature is reached.

Note:

If any invalid firing data were entered, the program will not start. An acoustic signal sounds and the display shows the invalid value flashing. Check entered firing data.

Button "R"

Press this button and enter the program number (01-99) to display the desired program.

After pressing this button, the program number display flashes for one minute. During this time, the two-digit program number must be entered via the numeric keypad.

Caution !

Program numbers lower than 10 must be entered with a leading zero (01-09).

Button Lift "↑"

Press the button "↑" to move the lift upwards.

Button Lift "↓"

Press the button "↓" to move the lift downwards.

At the end of the program, press the button "↓" in the lower lift position to start the individual rapid cooling.

Button "S"

Values in the display are saved to a program number.

After pressing this button, the program number display flashes for one minute. During this time, the two-digit program number must be entered via the numeric keypad.

Caution !

Program numbers lower than 10 must be entered with a leading zero (01-09).

Button "E"

Pressing this button opens the "Edit" function and the flashing cursor is displayed. The cursor flashes for one minute, allowing you to enter or alter data.

Use the cursor buttons to move the cursor in the display.

Button "C"

Press this button to:

- delete a wrong entry
- switch off rapid cooling
- confirm an error message
- abort the safety question when loading fixed programs
- abort a vacuum
- abort post-firing in the "View" function.

Cursor buttons ⬅️➡️↶️↷️

Press these buttons to move the cursor within the display.

Pressing and holding one of the buttons ⬅️ or ➡️ causes automatic movement of the cursor.

5.5 Menues

The menu provides all functions that cannot be executed directly. It includes the screen forms shown below, containing the menu name in the first line and a sub-menu selection or a parameter in the second line. Since only one sub-menu can be displayed at a time, use the buttons ➡️↶️ to bring the other sub-menus to the display. Press the button "E" to branch off into the selected sub-menu. Select a default setting with the buttons ⬅️➡️, or change a value with the button "C". Press the button ➡️ or ↶️ to accept the changes or press the button "M" repeatedly until the program is displayed again. The following section describes the function of the individual buttons for every remplate:

Settings:

Menu:
Configuration

Instrument parameters
Delete program memory
Load fixed programs

- ➡️ = next sub-menu
- ↶️ = previous sub-menu
- E = select sub-menu
- M = back to main menu

Configuration:
Contrast YY

Print	optional
Rapid cooling	Yes/No
Signal tone	Yes/No
Language	YYYY
Fahrenheit	Yes/No
Standby	Yes/No

- ➡️ = next parameter
- ↶️ = previous parameter
- ⬅️➡️ = select default or change contrast
- M = back to main menu

Instrument parameters:
Service Code ★★★★★

Op. hrs. muffle	YYYY
Op. hrs. pump	YYYY
Calibration	
offset	ZZ
Frequency	50/60 Hz

- ➡️ = next parameter
- ↶️ = previous parameter
- 0-9 and C = change values (★★★★★)
- C = delete values (YYYY, ZZ)
- M = back to main menu
- ⬅️ = increasing value (ZZ)
- ➡️ = decreasing value (ZZ)

Delete program memory

- ➡️ = next sub-menu
- ↶️ = previous sub-menu
- E = select sub-menu
- M = back to main menu

Delete all programs?
Jes = E No = C

- Safety question when deleting fixed programs.
- E = Delete and back to main menu
- C = back to main menu
- M = back to main menu

Load fixed programs

- ⇨ = next sub-menu
- ⇩ = previous sub-menu
- E = select sub-menu
- M = back to main menu

Overwrite no. 90–99?
Yes = E No = C

- Safety question when loading fixed programs.
- E = Load and back to main menu
 - C = back to main menu
 - M = back to main menu

When loading the fixed programs, the programs 90–99 will be overwritten.

If a currently loaded program is in that range, it will be overwritten as well.

The test program is displayed with the original data.

6 Create individual programs

6.1 Load program

As an example, load the free program 01 with the button “R” ⇨ 01.

01	0000	000	0000	c0020
	00	00	00,0	00,0 0

Note:

You have only 60 seconds to enter the individual firing parameters.

6.2 Enter preheating temperature

Press button “E”. The cursor automatically selects the preheating temperature. Use the numeric keypad 0–9 to enter e.g. 500.

01	0500	000	0000	c0020
	00	00	00,0	00,0 0

6.3 Enter heating rate

Press the cursor button “⇨”. The cursor selects the heating rate. Use the numeric keypad 0–9 to enter e.g. 80. The heating rate can be set to „max“ (140 °C/min) with the button „C“.

01	0500	080	0000	c0020
	00	00	00,0	00,0 0

6.4 Enter firing temperature

Press the cursor button “⇨”. The cursor selects the firing temperature. Use the numeric keypad 0–9 to enter e.g. 940.

01	0500	080	0940	c0020
	00	00	00,0	00,0 0

6.5 Enter drying time

Press the cursor button “⇨”. The cursor selects the drying time. Use the numeric keypad 0–9 to enter e.g. 6.

01	0500	080	0940	c0020
	06	00	00,0	00,0 0

6.6 Enter preheating time

Press the cursor button “⇨”. The cursor selects the preheating time. Use the numeric keypad 0–9 to enter e.g. 3.

01	0500	080	0940	c0020
	06	03	00,0	00,0 0

6.7 Enter vacuum level

Press the cursor button “⇨”. The cursor selects the vacuum level. Use the numeric keypad 0–9 to enter e.g. 50.

01	0500	080	0940	c0020
	06	03	50	00,0 00,0 0

6.8 Enter vacuum time

Press the cursor button “⇨”. The cursor selects the vacuum time. Use the numeric keypad 0–9 to enter e.g. 1.0 (fixed decimal point).

01	0500	080	0940	c0020
	06	03	50	01,0 00,0 0

Note:

The vacuum time must be at least 0.3 min shorter than the firing time.

6.9 Enter firing time

Press the **cursor button** “↔”. The cursor selects the firing time. Use the **numeric keypad 0–9** to enter e.g. 2.0 (fixed decimal point).

01	0500	080	0940	c0020
	06	03	50	01,0 02,0 0

Note:

The firing time must be at least 0.3 min longer than the vacuum time.

6.10 Save program

Save your program by pressing the **button “S”** and subsequently entering the program number, e.g. 01, with the numeric keypad.

6.11 Start program

Press the **button “start/stop”**. The green LED lights up and the total firing time is displayed in the in the lower left display.

01	0500	080	0940	c0020
0017	06	03	50	01,0 02,0 0

During firing, the following functions are available with the appropriate buttons:

- | | |
|--|--------------------------|
| a. Move lift (View function) | Lift buttons „↓↑“ |
| b. Abort program | Button start/stop |
| c. Abort vacuum or View function | Button “C” |
| d. Start editing (change parameters) | Button “E” |
| e. Show all theoretical parameter values | Button “M” |
| f. Delete errors | Button “C” |

6.12 Copy program

Copy your program by pressing the **button “S”** and subsequently entering the desired program number with the numeric keypad.

6.13 Special program

A saved program can be changed to a single special program by changing individual data and starting the program without pressing the **button S**.

7 Specialised functions

7.1 Standby

Standby can be preset in the sub-menu “Configuration”. In addition to the existing vacuum, this function allows maintaining a temperature of 120 °C in the firing chamber in “Off” state.

This function is also useful to keep moisture out of the firing chamber.

Display:

Standby	
120	120
(theoretical °C)	(actual °C)

7.2 Vacuum as special function

Vacuum can be preset in the sub-menu “Configuration”. This function allows maintaining a vacuum in the firing chamber in “Off” state.

This function is used to keep moisture out of the firing chamber.

Display:

8	OFF	8
(theoretical value 800 hPa)		(actual value 800 hPa)

7.3 Rapid cooling (manual)

If quick cooling is not preset in the configuration, it can be switched on manually after the program is completed and the firing chamber opened by pressing the lift button “↓”. In this case, the actual temperature must be greater than the theoretical preheating temperature. Press the button “C” to abort the quick cooling.

7.4 View function (soldering)

During firing without vacuum it is permissible to move the firing chamber downwards with the lift button “↓” and stop it with the lift button “↑”, in order inspect the workpiece. The firing time is stopped and the cursor appears on the firing temperature. Using the cursor buttons “↔”, and “↔”, the firing temperature can now be decreased or increased. Press the lift button “↑” to close the firing chamber again. The program continues. After the temperature has been set, post-firing is possible for up to three minutes. The post-firing time counts up from zero. Post-firing can be aborted at any time with the “start/stop” button.

7.5 Power failure backup

The Multimat C is equipped with a power failure backup. This backup activates as soon as there is a mains power failure in idle state or during a running firing program. If the downtime is shorter than 60 seconds, the program will continue to run. At the end of the program the following message appears: "F 10" (there was a power failure).

If the downtime is longer than 60 seconds, the program will be aborted and the message "F 10" is also displayed.

Check the firing result !

7.6 Testing the firing chamber temperature with the silver test

Calibration set "Silver test, manual", order no. D 03 532 803 (optional)

The temperature regulation accuracy is precisely set in the factory. If at any time checking the firing chamber temperature becomes necessary, you have the following possibility.

You require:

- 1 Dentsply sample carrier
- 1 piece of silver wire – diameter: 0.3, length: 37 mm.

Procedure:

- Thoroughly heat furnace for approx. 1 hour at 650°C.
- Place silver wire into Dentsply sample carrier.
- Enter the following data for the test program:

Preheating temperature	650 °C
Drying time	0 minutes
Preheating time	3 minutes
Vacuum level	0 hPa
Vacuum time	0 minutes
Firing time	1 minute
Cooling	0
Firing temperature	961 °C (melting point of the silver wire)
- Temperature increase rate max = full power
- Open firing chamber and place sample carrier with the silver wire centrally on the firing socket.
- Start test program.
- If after the test procedure the silver wire has started to melt along its upper third, i.e. a melting bead has formed on the surface of the wire, the temperature is exact to ± 2 °C.
- If the silver wire surface has not started to melt, then the firing chamber temperature is too low. In this case, repeat the test program – each time with a temperature increase of 3 °C – until the required melting effect on the surface of the silver wire is achieved.
- If the silver wire has melted altogether, then the firing chamber temperature is too high. In this case, repeat the test program – each time with a new length of silver wire and a temperature reduction of 3 °C – until the required melt effect on the surface of the silver wire is achieved.

If a deviation of the temperature from the theoretical value is detected, it can be corrected as follows:

1. Press button "Menu".
2. Press button \triangleleft to select sub-menu "Instrument parameters".
3. Confirm with button "E".
4. Press button \triangleleft to select sub-menu "Service Code".
5. Enter code no. 6070.
6. Press button \triangleleft to select sub-menu "Calibration offset".
7. Enter the desired correction value with the buttons \triangleleft (decrease value) and \triangleright (increase value).
8. Press sensor button "M" repeatedly until the program is displayed again.

The correction of the firing chamber temperature is completed.

Increase chamber temperature \triangleright increase calibration offset

Reduce chamber temperature \triangleleft reduce calibration offset

Example:

An increase of the calibration offset by the value 3 means an increase of the chamber temperature by approx. 3 °C.

Acoustic signals

Short signal: at each key press

Long signal: invalid entry

Triple signal: at end of program

7.7 Printer connection (Multimat C Art. no. D 03.022.230 V)

Printing the firing protocols

Note:

For printer connection you will require:

1 x null-modem cable V 24, serial
with one male and one female connection,

9 pin female

25 pin male

Configure the printer according to the following parameters:

4800 baud / 8 bit / 1 stop bit / no parity

Settings on the Multimat C furnace:

1. Press menu button “M”,
the screen shows:

Menu:
Configuration

2. Press button “E” and scroll up with the
cursor button \uparrow until “Print” is displayed,
the screen shows:

Configuration:
Print No

Now use the cursor button \rightarrow to change the
command Print “No” to Print “Yes”,
the screen shows:

Configuration:
Print Yes

3. Press menu button “M” twice,
the screen shows the last used program.

If the furnace was configured to “Print Yes”, it will send the firing
protocol to the printer interface after each firing.

If a printer with serial connection is connected and set up, the
sent firing protocol will be printed. If no printer is connected or
switched on, printing will be aborted.

8 Servicing and maintenance

8.1 Lift

The lift is used to vertically transport the workpiece.
Basically, it is a gear housing with rack and pinion. The unit is
driven by a synchronic geared motor.

If there is a power failure in the unit, the lift with the workpiece
can be pulled downwards by hand.

Maintenance note:

Apply a thin coating of nonresinous heat-resistant bearing
grease (local hardware store) on the rack, if required.

8.2 Firing plynth with firing plynth holder

The firing plynth is mounted vertically resilient over the firing
socket holder in the rack. It serves the purpose to insulate the
firing chamber on the bottom and simultaneously is used as
object carrier.

The firing plynth carrier with O-ring seals the firing chamber
during vacuum firing.

Maintenance note:

The seal of the firing plynth carrier with the O-ring must be kept
clean and must be inspected for possible damage regularly.

8.3 Vacuum pump

The Multimat C can be operated with any powerful vacuum
pump (power uptake max. 2.5 A), preferably with the Dentsply
vacuum pump. The setup is explained under section 4 “Setup
and first operation”. The pump connection tube must not be
longer than 2 meters.

Maintenance note:

Observe the operating manual of the pump for maintenance
purposes.

Important: Change oil every 3 months.

8.4 Replacing the firing muffle (trained electricians only)

Caution !

The heat insulation of the firing chamber in this
instrument consists of ceramic fibers. After extended
use of ceramic fibers at temperatures above 900 °C, silica
substances (cristobalite) may form. In certain cases, e.g. when
replacing the firing muffle, dust exposure may occur, which may
cause irritation of the skin, eyes, and respiratory organs.

Proceed as follows when replacing the firing muffle:

- All personnel must wear long sleeved clothing, headwear,
eye protection and gloves.
- Attach dust exhaustion to the dust source; if not possible,
equip personnel with dust mask FFP3.
- After completion of the tasks, possible adhering dust must
be first washed off the unprotected skin with cold water.
Afterwards wash with warm water and soap.
- Wash working clothes separately.

Warning !

**Disconnect the unit from main power before
opening the housing.**

Disassembly:

1. **Observe dust protection measures !**
2. Disconnect from main power.
3. Screw off cover and disconnect protective conductor.
4. Take off screws, cover and insulating disk.
5. Disconnect thermal element from both rear connecting bolts.
6. Disconnect heating wire ends of the muffle from both front
connecting bolts.
7. First remove the thermal element together with the bar, and
then the terminal insulation.
8. Remove muffler from the insulating insert; be careful if the
furnace is hot and touch only the ends of the heating wire.

Assembly:

Assemble with a new muffler in reverse order.

Caution!

Please wear cotton gloves when installing the new muffler, as to prevent sweat from the hands on the quartz tube.

Warning!

Check that the wire ends are not in contact with the metal mantle of the upper furnace section.

Shorten wire ends, if required. Reconnect protective conductor to the cover!

**8.5 Replacing the control element
(trained personnel only)**

1. **Disconnect from main power.**
 2. Unmount workpiece placement area.
 3. Unscrew both Phillips screws below the control element.
 4. Unscrew both Phillips screws above the control element.
 5. Swing control element outside.
 6. Disconnect vacuum tube and all plug connections.
 7. Disconnect thermal element.
 8. Disconnect protective conductor from the control unit.
- Installation of the new control unit is done in reverse order.

9 Faults and remedies

The following information shall assist you in identifying and evaluating simple errors and to fix them by appropriate remedial action.

A fault report consists of a two-digit number with a leading "F".

This report is on the very left in the second display row.

If several reports are generated simultaneously, the others are saved internally.

After confirming the displayed report with the button "C", the current line, or the next information message, is displayed again.

Warning!

Always disconnect furnace from mains power before opening the housing.

Fault:	Cause:	Remedy:
1. Display remains blank.	Mains fuse or transformer fuse faulty. Display faulty.	Disconnect from main power! Replace faulty fuses. If the display remains blank, contact service technician.
2. F 01 = Vacuum was not removed.	Ventilation valve faulty.	Contact service technician.
3. F 02 = Thermal el. ADC overflow.	Internal measurement error.	If it happens repeatedly, contact service technician.
4. F 03 = Chamber not fully closed.	Lift switch misaligned or faulty.	Contact service technician.
5. F 04 = Triac open.	Triac faulty.	Contact service technician.
6. F 05 = Muffle broken or Triac short circuit.	Muffle or Triac faulty.	Replace muffle or contact service technician.
7. F 06 = Control element too hot.	a. Missing cover plate over the control element. b. Extremely high standby temperature with open firing chamber.	Remove power plug. Allow furnace to cool for approx. 5 minutes. a. Mount cover plate. b. Keep furnace closed when not using it. Restart furnace normally after the cooling time.
8. F 07 = Theoretical vacuum was not reached.	Vacuum system leaks, suction performance of the pump too low or moisture in the firing chamber insulation.	Start program with vacuum. Observe vacuum display while vacuum pump is running. a. Display "9" = large leak, check sealing capability of the firing socket holder. b. Values 09–01 = Check vacuum system. c. Display between 99 hPa and theoretical value = Check vacuum system. If moisture has possibly entered the firing chamber insulation, heat to 1000 °C, hold for 10 min. and switch quick cooling on with opened firing chamber. d. Check vacuum pump according to manufacturer's instructions (change oil if required).

Fault:	Cause:	Remedy:
9. F 08 = Quick cooling is still on.	Idle temperature is not reached yet.	Wait until idle temperature is reached or abort manually with "C".
10. F 10 = There was a power failure.	Temporary power failure during firing.	Not possible. Please check the firing result.
11. F 15 = A vacuum is still present.	Vacuum system is not being ventilated.	Contact service technician.
12. F 16 = Thermal element broken.	Thermal element broken or interruption in thermal wiring.	Contact service technician.
13. F 17 = Thermal element poles reversed.	Plus and minus poles of thermal element or thermal wiring were interchanged.	Contact service technician.
14. F 18 = Thermal element short circuit.	Thermal element or thermal wiring contact each other.	Contact service technician.
15. F 19 = Overheated.	The actual temperature is 35 °C above the theoretical value.	Contact service technician.
16. F 20 = Program(s) faulty.	Memory error.	The control element deletes the faulty program automatically. The program must be entered again.

Caution!

Secure upper furnace part with transport safety prior to shipping the unit!

10 Preset programs

Firing table for Ceramco 3 metal ceramics (fixed programs assigned)

Pro-gram	Firing	Preheating temp. °C	Heating rate °C / min.	Firing temp. theoretical °C	Drying time min.	Preheating time min.	Vacuum level hPa	Vacuum time min.	Firing time min.	Cooling stage
92	Paste opaque	500	100	975	5	3	50	0.1	0.4	0
93	Powder opaque	650	70	970	3	3	50	0.1	0.4	0
94	Margin powder	650	70	965	5	5	50	0.1	0.4	0
95	Opaque dentine	650	45	930	5	5	50	0.1	1.0	0
96	Glaze w/o Glaze liquid	650	45	920	3	3	00	0	0.5	0
97	Glaze with Glaze liquid	650	55	925	3	3	00	0	0.5	0
98	Correcting powder	650	55	920	5	5	50	0.1	0.4	0

Firing table for Finesse metal ceramics (no fixed programs)

	Preheating temp. °C	Heating rate °C / min.	Firing temp. theoretical °C	Drying time min.	Preheating time min.	Vacuum level hPa	Vacuum time min.	Firing time min.
Powder Opaque	450	90	800	3	3	50	0.5	1
Paste Opaque	450	90	790	5	3	50	0.1	0.5
Correction Opaque	450	35	760	5	3	50	0.1	0.5
1. main firing	450	35	760	5	5	50	0.1	0.5
2. main firing	450	35	750	5	5	50	0.1	0.5
Correcting powder	450	55	730	5	5	50	0.1	0.5
Glaze	450	70	750	3	3	00	0	0.1

Firing temperatures are recommendations.

Make a test firing, if required, and adjust firing temperatures and times.

Set cooling stages according to alloy manufacturers' recommendations.

11 Warranty

Dentsply Furnaces

§ 1 Warranty scope

1. The manufacturer/warrantor provides the purchaser/warranty with a warranty that includes the functionality of the components mentioned in section 2 for the period of two years.
2. If such a component ceases to function within the warranty period, the purchaser/warranty is entitled to the required professional repair by replacing or repairing the component. The warranty does not include any claim to rescission of the sales contract or abatement (reduction of the purchasing price). The manufacturer repairs damage or faults of the unit that are based on manufacturing errors free of charge, provided they are reported immediately after identification.
3. The repairs included in the warranty also comprise testing, measuring and adjusting work, if necessary as a result or the remedy of a damage under warranty, but not maintenance work recommended by the manufacturer.

§ 2 Warranty extent

1. The warranty includes all assemblies and components of the unit.
2. In order to be entitled to warranty, the fault must be reported together with a copy of the invoice.
3. Claiming warranty does not extend the warranty period or start a new warranty period.
4. The warranty period for built in spare parts ends with the warranty period for the whole unit.

§ 3 Exclusion from warranty

1. The warranty expires, if persons who are not authorized by the manufacturer/warrantor carry out repairs or open the unit for purposes other than its intended use, or if the unit is fitted with spare parts that are not intended for this type of unit.
2. The warranty does not include transport damage, i.e. direct and sudden mechanical damage from outside.
3. The warranty expires if the unit is used improperly in ways not described in the operating manual or for other purposes.
4. The warranty does not cover claims beside those mentioned in section 1; in particular, it does not cover any damage arising outside of the unit, as long as liability is not required by law.

§ 4 Duties of the purchaser/warranty

1. The purchaser must have all maintenance, calibration and service work recommended in the operating manual/unit carried out at the depot/seller or in a workshop authorized by Dentsply.
2. The operating manual on the proper operation of the unit must be read prior to the first operation.
3. The special packing must be preserved carefully. If repairs become necessary, the unit must be packed according to the packing instructions.

§ 5 Transferring of claims and limitation of time

1. If the unit is sold while under warranty, the warranty claims are transferred to the new owner together with the ownership of the unit.
2. A warranty claim becomes prescribed 6 months after the damage occurred, and at the latest 6 months after the end of the duration of warranty.

§ 6 Consequences of a breach of duty

If the purchaser/warranty breaches one of the duties described in section 4 cl. 1 and 2, the manufacturer/warrantor is freed from their duty to provide warranty in this case.

§ 7 Legal liability

Legal liability remains unaffected by this warranty. This warranty exists beside the legally required warranty. If breaching the duty from section 4 cl. 3, the warranty expires, if the purchaser cannot prove that the breach of duty did not cause any damage.

§ 8 Applicable law

This warranty is subject to German law.
Exclusive court of jurisdiction is Hanau am Main.

Note:

Please also note the manufacturer's Terms and Conditions.

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we declare, that our product

Multimat C

Conforms with the requirements of the following Directive(s).

EU Council Directive 89/336/EEC, 73/23/EEC; changed to 93/68/EEC


The following harmonizing standards were applied:

**DIN EN 61010/03.94; DIN EN 50081-1: March 1993 in combination with
EN 55022: May 1995, EN 60555-2: 1997**

**DIN EN 50082-1: March 1993 in combination with IEC 801-2: 1984;
IEC 801-3: 1994, IEC 801-4: 1988**

Hanau, 10.07.2003

Hanau, 10.07.2003

Dr. Alexander Völcker
Manager Research and Development

Torsten Schwafert
Manager of Operations

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DeguDent GmbH