



MEDIRENT
your partner in clinical research

MORTARA ELI250

12 lead rest electrocardiograph



Technical data



Synopsis

Power Supply	100-240 V 50/60 Hz, 50 V Internally re-chargeable battery
Instrument Type	12-lead Electrocardiograph
Input Channels	Simultaneous acquisition of 12 leads
Standard Leads Acquired	I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6
Waveform Display	Backlit 1/4 VGA LCD, 3-channel lead groups or complete 8-channel presentation
Input Impedance	Meet or exceed the requirements of
Input Dynamic Range	ANSI/AAMI EC11
Electrode Offset Tolerance	
Common Mode Rejection	
ECG Storage	Up to 60 ECGs
Weight	5.1 kg including battery
Dimensions	39.4 x 43.2 x 10.2 cm

Price

Price will be based on the agreed minimum rental period and on specifications of required equipment.

Our standard minimum rental terms are 6, 12, 18, 24, 36, 48 and 60 months.

For a suitable tailor made quotation, please contact us at:

E: sales@medirent.com T: +31 71 581 581

Should you require different technical specifications, please do not hesitate to contact us.



MORTARA ELI250

12 lead rest electrocardiograph

Easy-to-use

The Mortara Eli 250 is an easy-to-use reliable 12-lead diagnostic electrocardiograph capable of viewing, acquiring, printing, and storing ECG data. Ideal for use in medical practices and small laboratories.

Backlit LCD Display

High-resolution display provides real time preview of 12-lead ECG and post-acquisition review of acquired ECG.

Full Keypad

Alphanumeric elastomer keypad also includes dedicated "one-touch" buttons for ECG acquisition, rhythm printing and ECG transmission.

Full-size, Z-fold Paper

Full-size paper (8.5" x 11" or A4 supported) output includes multiple Standard or Cabrera print options including 3+1, 3+3, 6-channel and 12-channel outputs.

Internal Storage

Internal storage for up to 60 digital ECG records. Optional expanded storage for up to 150 digital ECG records.

ECG quality control

This cardiograph provides extensive ECG quality-control features including AC noise filtering. The display allows preview of any three leads at a time to verify ECG quality before recording. The patient cables are well isolated to insure safety.



Contact information

Tieleman en Drosweg 17-19 / 2314 XW Leiden / The Netherlands
T +31 71 581 4 581 / F +31 71 581 4 582 / www.medirent.com